

WEST

Help

Logout

Interrupt

Main Menu

Search Form

Posting Counts

Show S Numbers

Edit S Numbers

Preferences

Cases

Search Results -

Term	Documents
(28 AND 24).USPT,PGPB.	0
(L28 AND L24).USPT,PGPB.	0

Database:

US Patents Full-Text Database
 US Pre-Grant Publication Full-Text Database
 JPO Abstracts Database
 EPO Abstracts Database
 Derwent World Patents Index
 IBM Technical Disclosure Bulletins

Search:

L32

Refine Search

Recall Text

Clear

Search History

DATE: Sunday, July 13, 2003 [Printable Copy](#) [Create Case](#)Set Name Query

side by side

Hit Count Set Name

result set

DB=USPT,PGPB; PLUR=YES; OP=ADJ

L32 L28 and L240 L32L31 (((435/91.1)!.CCLS.))2685 L31L30 (((435/270)!.CCLS.))172 L30

DB=USPT,PGPB,JPAB,EPAB,DWPI; PLUR=YES; OP=ADJ

L29 (((435/270)!.IPC.))0 L29

DB=USPT,PGPB; PLUR=YES; OP=ADJ

L28 (((536/25.4)!.CCLS.))606 L28

DB=USPT,PGPB,JPAB,EPAB,DWPI; PLUR=YES; OP=ADJ

L27 (((435/91.1)!.IPC.))0 L27L26 (((536/25.4)!.IPC.))0 L26L25 (((436/94)!.IPC.))0 L25L24 (silica gel) and L23286 L24

<u>L23</u>	chromatography and L21	310	<u>L23</u>
<u>L22</u>	(nucleic acid purification)and L20 and L19 and L18	2	<u>L22</u>
<u>L21</u>	L17 and L18 and L19 and L20	349	<u>L21</u>
<u>L20</u>	silica	400227	<u>L20</u>
<u>L19</u>	silicon dioxide	85653	<u>L19</u>
<u>L18</u>	xanthine	9837	<u>L18</u>
<u>L17</u>	nucleic acid	112449	<u>L17</u>

DB=USPT,PGPB; PLUR=YES; OP=ADJ

<u>L16</u>	L12 and l8	0	<u>L16</u>
<u>L15</u>	((435/91.1)!.CCLS.)	2685	<u>L15</u>
<u>L14</u>	((435/270)!.CCLS.)	172	<u>L14</u>

DB=USPT,PGPB,JPAB,EPAB,DWPI; PLUR=YES; OP=ADJ

<u>L13</u>	((435/270)!.IPC.)	0	<u>L13</u>
------------	---------------------	---	------------

DB=USPT,PGPB; PLUR=YES; OP=ADJ

<u>L12</u>	((536/25.4)!.CCLS.)	606	<u>L12</u>
------------	-----------------------	-----	------------

DB=USPT,PGPB,JPAB,EPAB,DWPI; PLUR=YES; OP=ADJ

<u>L11</u>	((435/91.1)!.IPC.)	0	<u>L11</u>
<u>L10</u>	((536/25.4)!.IPC.)	0	<u>L10</u>
<u>L9</u>	((436/94)!.IPC.)	0	<u>L9</u>
<u>L8</u>	(silica gel) and l7	286	<u>L8</u>
<u>L7</u>	chromatography and l5	310	<u>L7</u>
<u>L6</u>	(nucleic acid purification)and l4 and l3 and l2	2	<u>L6</u>
<u>L5</u>	l1 and l2 and l3 and L4	349	<u>L5</u>
<u>L4</u>	silica	400227	<u>L4</u>
<u>L3</u>	silicon dioxide	85653	<u>L3</u>
<u>L2</u>	xanthine	9837	<u>L2</u>
<u>L1</u>	nucleic acid	112449	<u>L1</u>

END OF SEARCH HISTORY

FILE 'BIOSIS, CAPLUS, BIOTECHNO' ENTERED AT 16:47:30 ON 13 JUL 2003

L1	437637 S SILICA
L2	1076 S L1 AND NUCLEIC ACID
L3	8 S L2 AND XANTHINE
L4	88477 S L1(W)GEL
L5	306 S L4 AND NUCLEIC ACID
L6	17 S L4 AND NUCLEIC ACID PURIFICATION
L7	0 S L6 AND XANTHINE
L8	4 S L5 AND XANTHINE
L9	2437 S CHROMATOGRAPHY AND XANTHINE
L10	110 S L9 AND SILICA
L11	74 S L9 AND SILICA GEL
L12	59 DUPLICATE REMOVE L11 (15 DUPLICATES REMOVED)
L13	4 S L10 AND NUCLEIC ACID
L14	4 S SILICON DIOXIDE AND L5
L15	4 S SILICON DIOXIDE AND NUCLEIC ACID PURIFICATION

L14 ANSWER 3 OF 4 CAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 2002:51620 CAPLUS

DOCUMENT NUMBER: 136:97266

TITLE: Isolating **nucleic acids** by
selective adsorption and desorption onto
silicon dioxide

INVENTOR(S): Weber, Martin; Singer, Thorsten; Cosaert, Sarah

PATENT ASSIGNEE(S): Qiagen G.m.b.H., Germany

SOURCE: PCT Int. Appl., 21 pp.

CODEN: PIXXD2

DOCUMENT TYPE: Patent

LANGUAGE: German

FAMILY ACC. NUM. COUNT: 1

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
WO 2002004620	A2	20020117	WO 2001-EP8066	20010712
WO 2002004620	A3	20020718		
W: JP, US				
RW: AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, TR				
DE 10033991	A1	20020124	DE 2000-10033991	20000712
EP 1299531	A2	20030409	EP 2001-971766	20010712
R: AT, BE, CH, DE, DK, ES, FR, GB, GR, IT, LI, LU, NL, SE, MC, PT, IE, FI, CY, TR				

PRIORITY APPLN. INFO.: DE 2000-10033991 A 20000712
WO 2001-EP8066 W 20010712

TI Isolating **nucleic acids** by selective adsorption and desorption onto **silicon dioxide**

AB The invention relates to a method for isolating **nucleic acids** from a soln., wherein the **nucleic acids** are adsorbed on a surface contg. SiO₂ in the presence of alkali halides and alc. The invention also relates to the use of a buffer soln. contg. alkali halides for isolating **nucleic acids** on a carrier contg. SiO₂, in addn. to a kit for implementing a method for isolating **nucleic acids** from a soln. The use of alkali metal halides avoids the use of hazardous chaotropic denaturants. Optimization expts. selecting appropriate salts and alcs. and ratios of salt to alc. are reported.

ST **nucleic acid** purifn silica sorbent; alkali halide alc
nucleic acid purifn; chloride isopropanol ethanol
nucleic acid purifn

IT Alcohols, uses
RL: MOA (Modifier or additive use); USES (Uses)
(C1-5; isolating **nucleic acids** by selective adsorption and desorption onto **silicon dioxide**)

IT Glass fibers, uses
Silica gel, uses
RL: DEV (Device component use); USES (Uses)
(as sorbent; isolating **nucleic acids** by selective adsorption and desorption onto **silicon dioxide**)

IT Plasmids
(isolating **nucleic acids** by selective adsorption and desorption onto **silicon dioxide**)

IT Alcohols, uses
Alkali metal halides, uses
RL: MOA (Modifier or additive use); USES (Uses)
(isolating **nucleic acids** by selective adsorption and desorption onto **silicon dioxide**)

IT **Nucleic acids**
RL: PUR (Purification or recovery); PREP (Preparation)
(isolating **nucleic acids** by selective adsorption and desorption onto **silicon dioxide**)

IT Synthetic fibers
 RL: DEV (Device component use); USES (Uses)
 (quartz, as sorbent; isolating **nucleic acids** by
 selective adsorption and desorption onto **silicon
 dioxide**)

IT 77-86-1, Tris (buffer) 1132-61-2; MOPS 7365-45-9, HEPES
 RL: MOA (Modifier or additive use); USES (Uses)
 (as buffer in elution medium; isolating **nucleic acids**
 by selective adsorption and desorption onto **silicon
 dioxide**)

IT 7631-86-9, **Silicon dioxide**, uses
 RL: DEV (Device component use); USES (Uses)
 (isolating **nucleic acids** by selective adsorption
 and desorption onto **silicon dioxide**)

IT 64-17-5, Ethanol, uses 67-63-0, Isopropanol, uses 7447-40-7, Potassium
 chloride, uses 7447-41-8, Lithium chloride, uses 7647-14-5, Sodium
 chloride, uses 8013-53-4
 RL: MOA (Modifier or additive use); USES (Uses)
 (isolating **nucleic acids** by selective adsorption
 and desorption onto **silicon dioxide**)

L14 ANSWER 4 OF 4 CAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 1988:403457 CAPLUS

DOCUMENT NUMBER: 109:3457

TITLE: Chromatographic process and resin preparation for the
 separation of **nucleic acids**

INVENTOR(S): Riesner, Detlev; Colpan, Metin

PATENT ASSIGNEE(S): Fed. Rep. Ger.

SOURCE: U.S., 8 pp. Cont.-in-part of U.S. Ser. No.560,931,
 abandoned.

CODEN: USXXAM

DOCUMENT TYPE: Patent

LANGUAGE: English

FAMILY ACC. NUM. COUNT: 2

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
US 4699717	A	19871013	US 1986-830708	19860214
DE 3211309	A1	19830929	DE 1982-3211309	19820326
DE 3211309	C2	19900816		

PRIORITY APPLN. INFO.: DE 1982-3211309 19820326
 US 1983-560931 19831125

TI Chromatographic process and resin preparation for the separation of
nucleic acids

AB A process for the chromatog. sepn. of **nucleic acids**
 using a chromatog. carrier material is described in which the surface of
 the carrier material is specially modified. **Silica gel**
 (particle diam. 10 .mu.m, pore size 4000 .ANG.) 50 g was heat and pressure
 activated, refluxed with .gamma.-glycidooxypropyltrimethoxysilane 100 mL
 in PhMe for 10 h under N, filtered, washed, and refluxed with
 N,N-diethylaminoethanol 100 mL and BF3/ether 1 mL for 12 h to give 51.5 g
 product. Potato spindle tuber viroid RNA from infected plants was
 chromatog. purified using the weak anion-exchange resin and a gradient
 elution of increasing KCl concn. in 5M urea, 30 mM K phosphate buffer, pH
 6.5. The purified **nucleic acid** was fully active in
 enzymic expts. (no data).

ST chromatog polymer stationary phase **nucleic acid**; plant
 virus RNA chromatog; aminosilane silica chromatog polynucleotide;
 silanized silica chromatog polynucleotide

IT Plasmid and Episome
 (DNA of, sepn. of, by chromatog. on aminosilane-modified **silica
 gel**)

IT Virus, plant

(RNA of, sepn. of, by chromatog. on aminosilane-modified **silica gel**)

IT Deoxyribonucleic acids
 RL: ANST (Analytical study)
 (sepn. of fragments of, by chromatog. on aminosilane-modified **silica gel**)

IT **Nucleic acids**
 Ribonucleic acids
 Ribonucleic acids, ribosomal
 Ribonucleic acids, transfer
 Ribonucleic acids, viral
 RL: PROC (Process)
 (sepn. of, by chromatog. on aminosilane-modified **silica gel**)

IT **Silica gel**, compounds
 RL: SPN (Synthetic preparation); PREP (Preparation)
 (aminosilylated, prepn. of, as chromatog. stationary phase for **nucleic acid sepn.**)

IT Nucleotides, polymers
 RL: PROC (Process)
 (oligo-, sepn. of, by chromatog. on aminosilane-modified **silica gel**)

IT Viroid
 (potato spindle tuber, RNA of, sepn. of, by chromatog. on aminosilane-modified **silica gel**)

IT Chromatography, column and liquid
 (stationary phases, aminosilane-modified **silica gel**, for **nucleic acid sepn.**)

IT 100-37-8DP, reaction products with glycidyloxypropyltrimethoxysilane and **silica gel** 2530-83-8DP, reaction products with diethylaminoethanol and **silica gel** 7631-86-9DP, **Silicon dioxide**, silanized and anion or cation exchanger-functionalized
 RL: SPN (Synthetic preparation); PREP (Preparation)
 (prepn. of, as chromatog. stationary phase for **nucleic acid sepn.**)

L15 ANSWER 2 OF 4 CAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 2000:366091 CAPLUS
DOCUMENT NUMBER: 133:13373
TITLE: Glass-coated particles for the purification of nucleic acids
INVENTOR(S): Harttig, Herbert; Riedling, Michael; Mennig, Martin; Schmidt, Helmut
PATENT ASSIGNEE(S): Institut fuer Neue Materialien Gem. G.m.b.H.; Germany; Roche Diagnostics G.m.b.H.
SOURCE: Ger. Offen., 8 pp.
CODEN: GWXXBX
DOCUMENT TYPE: Patent
LANGUAGE: German
FAMILY ACC. NUM. COUNT: 2
PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
DE 19854973	A1	20000531	DE 1998-19854973	19981130
WO 2000032762	A1	20000608	WO 1999-EP8996	19991123
W: AU, CA, JP, KR, NO, NZ, US				
RW: AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE				
EP 1144620	A1	20011017	EP 1999-959300	19991123
R: AT, BE, CH, DE, DK, ES, FR, GB, GR, IT, LI, LU, NL, SE, MC, PT, IE, FI				
JP 2002531084	T2	20020924	JP 2000-585393	19991123
NO 2001002620	A	20010727	NO 2001-2620	20010529
US 6545143	B1	20030408	US 2002-856737	20020109
US 2003125542	A1	20030703	US 2003-371375	20030220
PRIORITY APPLN. INFO.:		DE 1998-19854973 A 19981130		
		DE 1998-19855259 A 19981130		
		WO 1999-EP8996 W 19991123		
		US 2002-856737 A1 20020109		
ST	zinc aluminoborosilicate glass DNA RNA purifn; nucleic acid purifn zinc aluminoborosilicate glass			
IT	1303-86-2, Boron oxide, uses 1305-78-8, Calcium oxide, uses 1344-28-1, Aluminum oxide, uses 7631-86-9, Silicon dioxide , uses 12136-45-7, Potassium oxide, uses			
RL:	DEV (Device component use); USES (Uses) (glass-contg.; glass-coated particles for purifn. of nucleic acids)			

L15 ANSWER 3 OF 4 CAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 1997:172478 CAPLUS
DOCUMENT NUMBER: 126:168786
TITLE: Magnetic pigment
INVENTOR(S): Kleiber, Joerg; Walter, Thomas; Harttig, Herbert; Lesniak, Christoph; Mennig, Martin; Riedling, Michael; Schmidt, Helmut
PATENT ASSIGNEE(S): Boehringer Mannheim GmbH, Germany; Kleiber, Joerg; Walter, Thomas; Harttig, Herbert; Lesniak, Christoph; Mennig, Martin; Riedling, Michael; Schmidt, Helmut
SOURCE: PCT Int. Appl., 37 pp.
CODEN: PIXXD2
DOCUMENT TYPE: Patent
LANGUAGE: German
FAMILY ACC. NUM. COUNT: 2
PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
WO 9641811	A1	19961227	WO 1996-EP2459	19960606
W: AU, CA, CN, JP, KR, NO, NZ, US				
RW: AT, BE, CH, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE				

DE 19520398	A1	19961212	DE 1995-19520398	19950608
DE 19537985	A1	19970417	DE 1995-19537985	19951012
AU 9663007	A1	19970109	AU 1996-63007	19960606
AU 707115	B2	19990701		
EP 837871	A1	19980429	EP 1996-921935	19960606
EP 837871	B1	20030502		

R: AT, BE, CH, DE, DK, ES, FR, GB, IT, LI, NL, SE, IE, FI

JP 11509364	T2	19990817	JP 1997-502593	19960606
AT 239031	E	20030515	AT 1996-921935	19960606
NO 9705772	A	19980206	NO 1997-5772	19971208
US 6255477	B1	20010703	US 1998-952969	19980311

US-2002137920	A1	20020926	US 2001-756743	20010110
---------------	----	----------	----------------	----------

PRIORITY APPLN. INFO.:

DE 1995-19520398	A	19950608
DE 1995-19537985	A	19951012
WO 1996-EP2459	W	19960606
US 1998-952969	A3	19980311

ST biol material purifn magnetic pigment particle; **nucleic acid purifn** glass magnetic particle

IT 1303-86-2, Boron trioxide, uses 1309-38-2, Magnetite, uses 1314-23-4, Zirconium oxide, uses 1314-56-3, Phosphorus pentoxide, uses 1317-61-9, Iron oxide (Fe3O4), uses 1344-28-1, Aluminum oxide (Al2O3), uses 7631-86-9, **Silicon dioxide**, uses

RL: NUU (Other use, unclassified); USES (Uses)

(magnetic particles for purifn. of biol. materials)

L13 ANSWER 1 OF 4 CAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 2000:384565 CAPLUS

DOCUMENT NUMBER: 133:28236

TITLE: Methods and compositions for performing an array of chemical reactions on a support surface

INVENTOR(S): Zebala, John A.

PATENT ASSIGNEE(S): Syntrix Biochip, Inc., USA

SOURCE: PCT Int. Appl., 157 pp.

CODEN: PIXXD2

DOCUMENT TYPE: Patent

LANGUAGE: English

FAMILY-ACC.-NUM.-COUNT: 4

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
WO 2000033084	A2	20000608	WO 1999-US28021	19991123
WO 2000033084	A3	20000810		
W:	AE, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, CA, CH, CN, CR, CU, CZ, DE, DK, DM, EE, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MA, MD, MG, MK, MN, MW, MX, NO, NZ, PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TR, TT, TZ, UA, UG, US, UZ, VN, YU, ZA, ZW, AM, AZ, BY, KG, KZ, MD, RU, TJ, TM			
RW:	GH, GM, KE, LS, MW, SD, SL, SZ, TZ, UG, ZW, AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, BF, BJ, CF, CG, CI, CM, GA, GN, GW, ML, MR, NE, SN, TD, TG			
AU 2000018317	A5	20000619	AU 2000-18317	19991123
EP 1163374	A2	20011219	EP 1999-961813	19991123
R:	AT, BE, CH, DE, DK, ES, FR, GB, GR, IT, LI, LU, NL, SE, MC, PT, IE, SI, LT, LV, FI, RO			
JP 2002531470	T2	20020924	JP 2000-585669	19991123
PRIORITY APPLN. INFO.:			US 1998-110527P	P 19981201
			US 1999-326479	A 19990604
			WO 1999-US28021	W 19991123
ST	support array chem reaction photoresist; ligand array; DNA hybridization immobilized probe; ACE inhibitor screening enalaprilat analog solid phase synthesis; nucleic acid array			
IT	Nucleic acid hybridization (DNA-DNA; methods and compns. for performing arrays of chem. reactions on support surfaces using photoresists)			
IT	Adhesives Analysis Chromatography DNA sequence analysis Diagnosis Drug screening Electrophoresis Human immunodeficiency virus Indicators Mass spectrometry NMR spectroscopy Negative photoresists Nucleic acid hybridization PCR (polymerase chain reaction) Photoresists Positive photoresists Protein sequence analysis RNA sequence analysis Radiation Reactors Solvents Surface Synthesis			

(methods and compns. for performing arrays of chem. reactions on support surfaces using photoresists)

IT Probes (**nucleic acid**)
 RL: ARG (Analytical reagent use); ANST (Analytical study); USES (Uses)
 (methods and compns. for performing arrays of chem. reactions on support surfaces using photoresists)

IT Peptide **nucleic acids**
 RL: ARG (Analytical reagent use); DEV (Device component use); PEP (Physical, engineering or chemical process); RCT (Reactant); SPN (Synthetic preparation); ANST (Analytical study); PREP (Preparation); PROC (Process); RACT (Reactant or reagent); USES (Uses)
 (methods and compns. for performing arrays of chem. reactions on support surfaces using photoresists)

IT **Nucleic acids**
 Polynucleotides
 Proteins, general, reactions
 Reagents
 RL: ARG (Analytical reagent use); DEV (Device component use); RCT (Reactant); ANST (Analytical study); RACT (Reactant or reagent); USES (Uses)
 (methods and compns. for performing arrays of chem. reactions on support surfaces using photoresists)

IT Peptides, reactions
 RL: ARG (Analytical reagent use); DEV (Device component use); PEP (Physical, engineering or chemical process); RCT (Reactant); ANST (Analytical study); PROC (Process); RACT (Reactant or reagent); USES (Uses)
 (**nucleic acid** mimics; methods and compns. for performing arrays of chem. reactions on support surfaces using photoresists)

IT 51-20-7, 5-Bromouracil 51-21-8, 5-Fluorouracil 58-63-9, Inosine 65-71-4, Thymine 66-22-8, Uracil, uses 66-22-8D, Uracil, pseudo-, derivs., uses 68-94-0, Hypoxanthine 69-89-6, **Xanthine** 71-30-7, Cytosine 73-24-5, Adenine, uses 73-40-5, Guanine 141-90-2, Thiouracil 333-49-3, 2-Thiocytosine 443-72-1 504-07-4, Dihydrouracil 554-01-8, 5-Methylcytosine 578-76-7, 7-Methylguanine 591-28-6, 4-Thiouracil 636-26-0, 5-Methyl-2-thiouracil 696-07-1, 5-Iodouracil 938-85-2, 1-Methylguanine 1445-08-5, 2-Methyladenine 1445-15-4 1500-85-2, 7-Deazaadenine 1820-81-1, 5-Chlorouracil 1904-98-9, 2,6-Diaminopurine 2140-73-0, 1-Methylinosine 2365-40-4, N6-Isopentenyladenine 4776-08-3, 3-Methylcytosine 5142-22-3, 1-Methyladenine 6623-81-0, 5-Methoxyuracil 7355-55-7, 7-Deazaguanine 10030-78-1 14631-20-0 14886-75-0 20758-33-2 31458-37-4 72704-66-6 273752-46-8 273752-47-9 273752-48-0 273752-50-4 273752-52-6
 RL: DEV (Device component use); PRP (Properties); USES (Uses)
 (array of nucleobase polymers contg.; methods and compns. for performing arrays of chem. reactions on support surfaces using photoresists)

IT 7631-86-9, **Silica**, uses
 RL: DEV (Device component use); USES (Uses)
 (methods and compns. for performing arrays of chem. reactions on support surfaces using photoresists)

L13 ANSWER 2 OF 4 CAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 1995:507908 CAPLUS

DOCUMENT NUMBER: 122:265933

TITLE: Preparation of pyranose nucleoside derivatives as antiviral and antitumor agents

INVENTOR(S): Waga, Toshiaki; Meguro, Hiromu; Oorui, Hiroshi

PATENT ASSIGNEE(S): Asahi Breweries Ltd, Japan

SOURCE: Jpn. Kokai Tokkyo Koho, 10 pp.

CODEN: JKXXAF

DOCUMENT TYPE: Patent

LANGUAGE: Japanese
FAMILY ACC. NUM. COUNT: 1
PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
JP 06263793	A2	19940920	JP 1993-139791	19930304
PRIORITY APPLN. INFO.:			JP 1993-139791	19930304

OTHER SOURCE(S): MARPAT 122:265933

AB The title compds. (I; B = adenine, guanine, thymine, uracil, cytosine, hypoxanthine, **xanthine**, 5-methylcytosine, 4-ethoxy-5-methyl-2-oxopyrimidine, 4-isopropoxy-5-methyl-2-oxopyrimidine, 5-methyl-2-oxopyrimidine; R1, R2 = H, OH; or R1R2 = bond; R3 = Q wherein n = 0,1,3; R4 = H, lower alkoxy) or pharmacol. acceptable esters, ethers, or salts thereof are prepd. as antiviral and antitumor agents, particularly potential anti-HIV agents (no data), are prepd. Thus, 2.0 g adenine and 2.0 g K2CO3 were suspended in 100 mL DMF and after stirring at 80.degree. for 1 h, 2.0 g 18-crown-6 ether and Me 2,3-anhydro-4,6-O-benzylidene-.alpha.-D-allopyranoside were added followed by stirring the resulting mixt. at 120.degree. for 16 h to give, after **silica gel chromatog.**, 91% adenylaltropyranoside deriv. (II; RR = CHPh).

IT 3150-15-0

RL: RCT (Reactant); RACT (Reactant or reagent)
(condensation with **nucleic acid** bases in prepn. of
pyranose nucleoside derivs. as antiviral and antitumor agents)

L13 ANSWER 3 OF 4 CAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 1984:628477 CAPLUS

DOCUMENT NUMBER: 101:228477

TITLE: Separation of nucleobases on polar amino cyano
high-performance liquid **chromatography**
columns

AUTHOR(S): Joshua, Henry; Goetz, Michael

CORPORATE SOURCE: Merck Sharp and Dohme Res. Lab., Rahway, NJ, 07065,
USA

SOURCE: Journal of Chromatography (1984), 303(1), 185-9
CODEN: JOCRAM; ISSN: 0021-9673

DOCUMENT TYPE: Journal

LANGUAGE: English

TI Separation of nucleobases on polar amino cyano high-performance liquid
chromatography columns

AB The use of a polar amino cyano (PAC) HPLC column in conjunction with buffered water-acetonitrile (8:92) eluents affords a viable alternative to the std. reversed-phase, **silica gel** and ion-exchange methods for the **chromatog.** sepn. of nucleobases. A Whatman PAC column (Partisil PXS 5/25 PAC, 25 cm .times. 4.6 mm, 5 .mu.m particle size) was used equipped with an Upchurch C-130 precolumn packed with Whatman Co-Pell PAC 30-38 .mu.m particles. Seven nucleobases (adenine [73-24-5], cytosine [71-30-7], guanine [73-40-5], hypoxanthine [68-94-0], thymine [65-71-4], uracil [66-22-8], and **xanthine** [69-89-6]) were successfully sepd. Changes in the eluent pH values were found to affect the selectivity and capacity factors for the nucleobases. Thus effects were esp. pronounced with **xanthine**. The method is useful for the sepn. of nucleobases from fermn. broth exts.

ST **nucleic acid** base sepn polar HPLC; amino cyano
chromatog adenine cytosine guanine

IT Fermentation

(**nucleic acid** bases, sepn. after, by polar
aminocyano high-performance liq. **chromatog.**)

IT **Nucleic acids**

RL: PROC (Process)

(bases, sepn. of, from fermn. broth by polar amino cyano
high-performance liq. **chromatog.**)

IT **Chromatography**, column and liquid

(high-performance, **nucleic acid** bases sepn. from
fermn. broth by, on polar amino cyano stationary phase)
IT 65-71-4 66-22-8, analysis 68-94-0 69-89-6 71-30-7 73-24-5,
analysis 73-40-5
RL: PROC (Process)
(sepn. of, from fermn. broth by polar amino cyano high-performance liq.
chromatog.)

L13 ANSWER 4 OF 4 CAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 1968:48556 CAPLUS

DOCUMENT NUMBER: 68:48556

TITLE: Mode of action of the chemosterilants,
2-imidazolidinone and 4-imidazolin-2-one, in the
housefly and in the large milkweed bug

AUTHOR(S): Schaefer, Charles Herbert

CORPORATE SOURCE: Shell Develop. Co., Modesto, CA, USA

SOURCE: Life Sciences (1967), 6(24), 2677-83

CODEN: LIFSAK; ISSN: 0024-3205

DOCUMENT TYPE: Journal

LANGUAGE: English

AB An anal. method was developed for detecting 2-imidazolidinone (I) and
4-imidazolin-2-one (II) in insect tissues and feces after injection of 5
.mu.g. I or II into female houseflies or large milkweed bugs, *Oncopeltus*
fasciatus. Insects were homogenized and the homogenate or feces were
extd. with MeOH. After evapn., the residue was added to a H₂O-CHCl₃
mixt., centrifuged, the aq. phase extd. with CHCl₃, sepd., and evapd. onto
silica gel. Exts. were subjected to column **chromatog.**
and thin-layer **chromatog.** Excretion plus metabolism eliminated
84% of a 5-.mu.g. dose of I and 76% of II within 24 hrs. after injection
into houseflies; these 2 compds. are temporary sterilants in this species.
The nature of the metabolites of either compd. was not detd. Large
milkweed bugs were apparently unable to detoxify either compd. and there
was no trace of either in the feces at 48 hrs. after injection; both I and
II produce permanent sterility in this species. An attempt was made to
det. the mode of action of II in houseflies by feeding the lowest level
that inhibited reproduction in the presence of potential reversers; none
of the natural biochem. intermediates tested (vitamins, glycine,
histamine, inosinic acid, orotic acid, cytidine, adenine, cytidylic acid,
adenylic acid, deoxyadenosine, **xanthine**, RNA, oleic acid,
.beta.-sitosterol, and cholesterol) had any effect on the sterilant
activity of II. No synergism was apparent when either compd. was fed in
the presence of 0.2% sesamex or when II was given in the presence of 1%
H₃BO₃. II may inhibit the formation of complex mols. such as proteins or
nucleic acids rather than that of simple mols.; the mode
of action of II is complex and may also involve interference with
endocrine regulation.

FILE 'HOME' ENTERED AT 16:47:21 ON 13 JUL 2003

=> FIL BIOSIS, CAPLUS, BIOTECHNO
COST IN U.S. DOLLARS

SINCE FILE ENTRY	TOTAL SESSION
0.21	0.21

FULL ESTIMATED COST

FILE 'BIOSIS' ENTERED AT 16:47:30 ON 13 JUL 2003

COPYRIGHT (C) 2003 BIOLOGICAL ABSTRACTS INC.(R)

FILE 'CAPLUS' ENTERED AT 16:47:30 ON 13 JUL 2003

USE IS SUBJECT TO THE TERMS OF YOUR STN CUSTOMER AGREEMENT.

PLEASE SEE "HELP USAGETERMS" FOR DETAILS.

COPYRIGHT (C) 2003 AMERICAN CHEMICAL SOCIETY (ACS)

FILE 'BIOTECHNO' ENTERED AT 16:47:30 ON 13 JUL 2003

COPYRIGHT (C) 2003 Elsevier Science B.V., Amsterdam. All rights reserved.

=> s silica

L1 437637 SILICA

=> s l1 and nucleic acid

2 FILES SEARCHED...

L2 1076 L1 AND NUCLEIC ACID

=> s l2 and xanthine

L3 8 L2 AND XANTHINE

=> d ibib kwic 1-8 l3

L3 ANSWER 1 OF 8 CAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 2002:90055 CAPLUS

DOCUMENT NUMBER: 136:131252

TITLE: Cationic materials and methods for covalent bonding
nucleic acids to high purity
silica surfaces

INVENTOR(S): Lyles, Mark B.

PATENT ASSIGNEE(S): USA

SOURCE: PCT Int. Appl., 9 pp.

CODEN: PIXXD2

DOCUMENT TYPE: Patent

LANGUAGE: English

FAMILY ACC. NUM. COUNT: 1

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
WO 2002008237	A2	20020131	WO 2001-US23079	20010720
WO 2002008237	A3	20021107		
W: AE, AG, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, BZ, CA, CH, CN, CO, CR, CU, CZ, DE, DK, DM, DZ, EC, EE, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MA, MD, MG, MK, MN, MW, MX, MZ, NO, NZ, PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TR, TT, TZ, UA, UG, UZ, VN, YU, ZA, ZW, AM, AZ, BY, KG, KZ, MD, RU, TJ, TM				
RW: GH, GM, KE, LS, MW, MZ, SD, SL, SZ, TZ, UG, ZW, AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, TR, BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, ML, MR, NE, SN, TD, TG				
AU 2001076023	A5	20020205	AU 2001-76023	20010720
US 2002103350	A1	20020801	US 2001-910697	20010720
EP 1305328	A2	20030502	EP 2001-953590	20010720
R: AT, BE, CH, DE, DK, ES, FR, GB, GR, IT, LI, LU, NL, SE, MC, PT, IE, SI, LT, LV, FI, RO, MK, CY, AL, TR				

Welcome to STN International! Enter x:x

LOGINID:SSSPTA1623HXO

PASSWORD:

TERMINAL (ENTER 1, 2, 3, OR ?):2

-----* * * * * Welcome to STN International * * * * *-----

NEWS	1	Web Page URLs for STN Seminar Schedule - N. America
NEWS	2	"Ask CAS" for self-help around the clock
NEWS	3	Feb 24 PCTGEN now available on STN
NEWS	4	Feb 24 TEMA now available on STN
NEWS	5	Feb 26 NTIS now allows simultaneous left and right truncation
NEWS	6	Feb 26 PCTFULL now contains images
NEWS	7	Mar 04 SDI PACKAGE for monthly delivery of multifile SDI results
NEWS	8	Mar 24 PATDPAFULL now available on STN
NEWS	9	Mar 24 Additional information for trade-named substances without structures available in REGISTRY
NEWS	10	Apr 11 Display formats in DGENE enhanced
NEWS	11	Apr 14 MEDLINE Reload
NEWS	12	Apr 17 Polymer searching in REGISTRY enhanced
NEWS	13	Jun 13 Indexing from 1947 to 1956 added to records in CA/CAPLUS
NEWS	14	Apr 21 New current-awareness alert (SDI) frequency in WPIDS/WPINDEX/WPIX
NEWS	15	Apr 28 RDISCLOSURE now available on STN
NEWS	16	May 05 Pharmacokinetic information and systematic chemical names added to PHAR
NEWS	17	May 15 MEDLINE file segment of TOXCENTER reloaded
NEWS	18	May 15 Supporter information for ENCOMPPAT and ENCOMPLIT updated
NEWS	19	May 19 Simultaneous left and right truncation added to WSCA
NEWS	20	May 19 RAPRA enhanced with new search field, simultaneous left and right truncation
NEWS	21	Jun 06 Simultaneous left and right truncation added to CBNB
NEWS	22	Jun 06 PASCAL enhanced with additional data
NEWS	23	Jun 20 2003 edition of the FSTA Thesaurus is now available
NEWS	24	Jun 25 HSDB has been reloaded
NEWS EXPRESS		April 4 CURRENT WINDOWS VERSION IS V6.01a, CURRENT MACINTOSH VERSION IS V6.0b(ENG) AND V6.0Jb(JP), AND CURRENT DISCOVER FILE IS DATED 01 APRIL 2003
NEWS HOURS		STN Operating Hours Plus Help Desk Availability
NEWS INTER		General Internet Information
NEWS LOGIN		Welcome Banner and News Items
NEWS PHONE		Direct Dial and Telecommunication Network Access to STN
NEWS WWW		CAS World Wide Web Site (general information)

Enter NEWS followed by the item number or name to see news on that specific topic.

All use of STN is subject to the provisions of the STN Customer agreement. Please note that this agreement limits use to scientific research. Use for software development or design or implementation of commercial gateways or other similar uses is prohibited and may result in loss of user privileges and other penalties.

* * * * * STN Columbus * * * * *

PRIORITY APPLN. INFO.:

US 2000-220096P P 20000721

US 2001-910697 A 20010720

WO 2001-US23079 W 20010720

- TI Cationic materials and methods for covalent bonding **nucleic acids** to high purity **silica** surfaces
- AB Surfaces contg. high purity **silica** (silicon dioxide) exhibit high loading potential for **nucleic acids**. Formulations contg. **nucleic acids** and materials which mask the electrostatic interactions between the **nucleic acids** and surfaces are disclosed. By masking the phosphate charges of the **nucleic acids**, undesired interactions may be minimized or eliminated, thereby allowing the covalent bonding of the **nucleic acids** to the surface to proceed. The use of such formulations addnl. minimizes nonspecific binding of the **nucleic acids** to the surface. Examples of materials to be included in such formulations include cations, **xanthines**, hexoses, purines, arginine, lysine, polyarginine, polylysine, and quaternary ammonium salts.
- ST **nucleic acid** covalent bond **silica** surface
electrostatic interaction cation
- IT Electrostatic force
(attractive, minimization of, by phosphate group masking; cationic materials and methods for covalent bonding **nucleic acids** to high purity **silica** surfaces)
- IT Spheres
(beads; cationic materials and methods for covalent bonding **nucleic acids** to high purity **silica** surfaces)
- IT Cations
Immobilization, molecular
(cationic materials and methods for covalent bonding **nucleic acids** to high purity **silica** surfaces)
- IT Hexoses
Quaternary ammonium compounds, uses
RL: NUU (Other use, unclassified); USES (Uses)
(cationic materials and methods for covalent bonding **nucleic acids** to high purity **silica** surfaces)
- IT Glass beads
RL: NUU (Other use, unclassified); RCT (Reactant); RACT (Reactant or reagent); USES (Uses)
(cationic materials and methods for covalent bonding **nucleic acids** to high purity **silica** surfaces)
- IT Bond
(covalent; cationic materials and methods for covalent bonding **nucleic acids** to high purity **silica** surfaces)
- IT Attractive force
(electrostatic, minimization of, by phosphate group masking; cationic materials and methods for covalent bonding **nucleic acids** to high purity **silica** surfaces)
- IT DNA
Nucleic acids
RL: PEP (Physical, engineering or chemical process); RCT (Reactant); PROC (Process); RACT (Reactant or reagent)
(immobilization of; cationic materials and methods for covalent bonding **nucleic acids** to high purity **silica** surfaces)
- IT Oligonucleotides
RL: SPN (Synthetic preparation); PREP (Preparation)
(immobilized; cationic materials and methods for covalent bonding **nucleic acids** to high purity **silica** surfaces)
- IT Phosphate group
(masking of phosphate groups to reduce nonspecific binding; cationic

materials and methods for covalent bonding **nucleic acids** to high purity **silica** surfaces)

IT Molecular association

(nonspecific, minimizing of; cationic materials and methods for covalent bonding **nucleic acids** to high purity **silica** surfaces)

IT Synthetic fibers

RL: NUU (Other use, unclassified); RCT (Reactant); RACT (Reactant or reagent); USES (Uses)

(**silica**; cationic materials and methods for covalent bonding **nucleic acids** to high purity **silica** surfaces)

IT 50-44-2, 6-Thiohypoxanthine 50-99-7, Glucose, uses 56-37-1, Benzyltriethyl ammonium chloride 56-87-1, Lysine, uses 56-93-9, Benzyltrimethyl ammonium chloride 57-48-7, Fructose, uses 58-08-2, Caffeine, uses 58-63-9, Inosine 59-23-4, Galactose, uses 68-94-0, Hypoxanthine 69-89-6, **Xanthine** 74-79-3, L-Arginine, uses 75-57-0, Tetramethyl ammonium chloride 83-67-0, 3,7-Dimethylxanthine 87-79-6, Sorbose 104-74-5, N-Lauryl pyridinium chloride 120-73-0, Purine 146-80-5, Xanthosine 519-32-4, 1,3,9-Trimethylxanthine 552-62-5, 7-Methylxanthine 574-25-4 611-59-6, 1,7-Dimethylxanthine 628-13-7D, Pyridinium chloride, N-alkyl derivs. 652-37-9 890-38-0, 2'-Deoxyinosine 1076-22-8, 3-Methylxanthine 1198-33-0, 9-Methylxanthine 1643-19-2, Tetrabutyl ammonium bromide 2002-59-7, 6-Thioxanthine 2036-13-7, 6-Purinecarbonitrile 3458-28-4, Mannose 5137-55-3, Trioctylmethyl ammonium chloride 5270-30-4 5437-25-2, 2,6-Dithiopurine 5438-71-1, 8-(3-Carboxypropyl)-1,3-dimethylxanthine 5987-68-8, Altrose 6038-51-3, Allose 6136-37-4, 1-Methylxanthine 6739-64-6, Nicotinamide hypoxanthine dinucleotide phosphate 14114-46-6, 3,7-Dimethyl-1-propargylxanthine 15837-08-8, 3,9-Dimethylxanthine 17338-96-4, 8-Methylxanthine 17598-81-1, Tagatose 23616-79-7, Benzyltributyl ammonium chloride 24937-47-1, Polyarginine 25104-18-1, Polylysine 25212-18-4, Polyarginine 28822-58-4, 3-Isobutyl-1-methylxanthine 30077-17-9, Talose 31542-51-5 31542-63-9, 1,3-Dipropyl-7-methylxanthine 31617-39-7, 1,3-Diethyl-7-methylxanthine 32503-27-8, Tetrabutyl ammonium hydrogensulfate 33073-01-7, 1,9-Dimethylxanthine 35873-49-5, 8-Cyclopentyl-1,3-dimethylxanthine 38000-06-5, Polylysine 41078-02-8, 3-Propylxanthine 55242-55-2, 3-Methyl-1-(5-oxohexyl)-7-propylxanthine 59840-67-4 70332-31-9 75922-48-4, 1,3-Diethyl-8-phenylxanthine 78033-08-6, 8-Methoxymethyl-3-isobutyl-1-methylxanthine 89073-57-4, 1,3-Dipropyl-8-p-sulfophenylxanthine 91725-06-3 96654-24-9 102146-07-6, 8-Cyclopentyl-1,3-dipropylxanthine 103258-00-0 135462-23-6 149981-23-7 149981-25-9 194802-32-9 195522-91-9 197456-29-4 199190-66-4 392687-18-2 392687-19-3 392687-20-6 392687-21-7

RL: NUU (Other use, unclassified); USES (Uses)

(cationic materials and methods for covalent bonding **nucleic acids** to high purity **silica** surfaces)

IT 7631-86-9, **Silica**, reactions

RL: NUU (Other use, unclassified); RCT (Reactant); RACT (Reactant or reagent); USES (Uses)

(cationic materials and methods for covalent bonding **nucleic acids** to high purity **silica** surfaces)

L3 ANSWER 2 OF 8 CAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 2001:284148 CAPLUS

DOCUMENT NUMBER: 134:306118

TITLE: Template-dependent ligation with PNA-DNA chimeric probes

INVENTOR(S): Egholm, Michael; Chen, Caifu

PATENT ASSIGNEE(S): PE Corporation, USA

SOURCE: PCT Int. Appl., 66 pp.

CODEN: PIXXD2

DOCUMENT TYPE: Patent
 LANGUAGE: English
 FAMILY ACC. NUM. COUNT: 1
 PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
WO 2001027326	A2	20010419	WO 2000-US27730	20001006
WO 2001027326	A3	20020510		
W: AE, AG, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, BZ, CA, CH, CN, CR, CU, CZ, DE, DK, DM, DZ, EE, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MA, MD, MG, MK, MN, MW, MX, MZ, NO, NZ, PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TR, TT, TZ, UA, UG, UZ, VN, YU, ZA, ZW, AM, AZ, BY, KG, KZ, MD, RU, TJ, TM RW: GH, GM, KE, LS, MW, MZ, SD, SL, SZ, TZ, UG, ZW, AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, BF, BJ, CF, CG, CI, CM, GA, GN, GW, ML, MR, NE, SN, TD, TG US 6297016 B1 20011002 US 1999-416003 19991008 EP 1220953 A2 20020710 EP 2000-968853 20001006 R: AT, BE, CH, DE, DK, ES, FR, GB, GR, IT, LI, LU, NL, SE, MC, PT, IE, SI, LT, LV, FI, RO, MK, CY, AL JP 2003511059 T2 20030325 JP 2001-529456 20001006 US 6469151 B1 20021022 US 2001-881557 20010614 US 2002177133 A1 20021128				

PRIORITY APPLN. INFO.: US 1999-416003 A 19991008
 WO 2000-US27730 W 20001006

AB The invention provides methods, kits, and compns. for ligation of peptide-
nucleic acid (PNA)-DNA chimeric probes and
 oligonucleotides when they are hybridized adjacently to template
nucleic acids using ligases and ligation reagents. The
 invention is based in part on the discovery that a ligase enzyme can
 ligate a PNA-DNA chimeric probe and a second probe under a broad range of
 exptl. conditions and variables. Structural requirements of the chimeras
 for ligation include 5 to 15 contiguous PNA monomer units, 2 or more
 contiguous nucleotides, and a 3' hydroxyl or 5' hydroxyl terminus. The
 chimera and/or oligonucleotide may be labeled with fluorescent dyes or
 other labels. The methods include, for example, oligonucleotide-ligation
 assays (OLA) and single nucleotide polymorphism detection.

ST PNA DNA chimeric probe ligation oligonucleotide; peptide **nucleic**
acid DNA chimera ligation oligonucleotide

IT Genetic methods
 (OLA (oligonucleotide ligation assay); template-dependent ligation with
 PNA (peptide-**nucleic acid**)-DNA chimeric probes)

IT Cyanine dyes
 (fluorescent dye label; template-dependent ligation with PNA (peptide-
nucleic acid)-DNA chimeric probes)

IT Genetic polymorphism
 (single nucleotide, detection of; template-dependent ligation with PNA
 (peptide-**nucleic acid**)-DNA chimeric probes)

IT Glass, uses
 Polyamides, uses
 Silica gel, uses

RL: TEM (Technical or engineered material use); USES (Uses)
 (solid support; template-dependent ligation with PNA (peptide-
nucleic acid)-DNA chimeric probes)

IT Test kits
 (template-dependent ligation with PNA (peptide-**nucleic**
acid)-DNA chimeric probes)

IT DNA
 Oligodeoxyribonucleotides
 Peptide **nucleic acids**
 Probes (**nucleic acid**)
 RNA

RL: BPR (Biological process); BSU (Biological study, unclassified); BUU (Biological use, unclassified); BIOL (Biological study); PROC (Process); USES (Uses)

(template-dependent ligation with PNA (peptide-**nucleic acid**)-DNA chimeric probes)

IT 76823-03-5, FAM

RL: BUU (Biological use, unclassified); BIOL (Biological study); USES (Uses)

(FAM, fluorescent dye label; template-dependent ligation with PNA (peptide-**nucleic acid**)-DNA chimeric probes)

IT 155911-16-3, HEX

RL: BUU (Biological use, unclassified); BIOL (Biological study); USES (Uses)

(HEX, fluorescent dye label; template-dependent ligation with PNA (peptide-**nucleic acid**)-DNA chimeric probes)

IT 51-28-5, 2,4-Dinitrophenol, biological studies 57-88-5, Cholesterol, biological studies 58-85-5, Biotin 1672-46-4, Digoxigenin 2321-07-5, Fluorescein

RL: BUU (Biological use, unclassified); BIOL (Biological study); USES (Uses)

(affinity ligand; template-dependent ligation with PNA (peptide-**nucleic acid**)-DNA chimeric probes)

IT 69-89-6D, **Xanthine**, arom.-substituted 596-12-3 82855-40-1, JOE 120718-39-0, ROX 120718-52-7, TAMRA 192230-82-3, TET 278175-13-6

RL: BUU (Biological use, unclassified); BIOL (Biological study); USES (Uses)

(fluorescent dye label; template-dependent ligation with PNA (peptide-**nucleic acid**)-DNA chimeric probes)

IT 569-64-2, Malachite Green 6268-49-1 56512-49-3, DABSYL chloride 245106-84-7, NTB 335277-36-6, d-TAMRA

RL: BUU (Biological use, unclassified); BIOL (Biological study); USES (Uses)

(fluorescent quencher; template-dependent ligation with PNA (peptide-**nucleic acid**)-DNA chimeric probes)

IT 1438-30-8, Netropsin 23491-45-4, Hoechst 33258 39389-47-4, Distamycin 114309-58-9

RL: BUU (Biological use, unclassified); BIOL (Biological study); USES (Uses)

(minor groove binder; template-dependent ligation with PNA (peptide-**nucleic acid**)-DNA chimeric probes)

IT 868-77-9 7631-86-9, **Silica**, uses 9002-88-4, Polyethylene 9003-01-4, Polyacrylic acid 9003-05-8, Polyacrylamide

RL: TEM (Technical or engineered material use); USES (Uses)

(solid support; template-dependent ligation with PNA (peptide-**nucleic acid**)-DNA chimeric probes)

IT 65-71-4D, Thymine, PNA-DNA chimera contg. 66-22-8D, Uracil, PNA-DNA chimera contg., biological studies 68-94-0D, Hypoxanthine, PNA-DNA chimera contg. 71-30-7D, Cytosine, PNA-DNA chimera contg. 73-24-5D, Adenine, PNA-DNA chimera contg. 73-40-5D, Guanine, PNA-DNA chimera contg. 135-67-1D, Phenoxazine, PNA-DNA chimera contg. 489-59-8D, Isocytidine, PNA-DNA chimera contg. 1450-85-7D, 2-Thiopyrimidine, PNA-DNA chimera contg. 1500-85-2D, 7-Deazaadenine, PNA-DNA chimera contg. 1818-71-9D, Isoguanosine, PNA-DNA chimera contg. 1904-98-9D, 2,6-Diaminopurine, PNA-DNA chimera contg. 4562-27-0D, 4(3H)-Pyrimidone, PNA-DNA chimera contg. 7355-55-7D, 7-Deazaguanine, PNA-DNA chimera contg. 13230-97-2D, 8-Oxopurine, PNA-DNA chimera contg. 24123-14-6D, N-2-Aminoethylglycine, PNA-DNA chimera contg. 57100-18-2D, Pseudoisocytidine, PNA-DNA chimera contg. 335084-02-1D, PNA-DNA chimera contg.

RL: BPR (Biological process); BSU (Biological study, unclassified); BUU (Biological use, unclassified); BIOL (Biological study); PROC (Process); USES (Uses)

(template-dependent ligation with PNA (peptide-**nucleic acid**)-DNA chimeric probes)

acid)-DNA chimeric probes)
IT 9015-85-4, DNA ligase 37211-65-7, Polynucleotide kinase 37353-39-2,
RNA ligase
RL: CAT (Catalyst use); USES (Uses)
(template-dependent ligation with PNA (peptide-nucleic
acid)-DNA chimeric probes)

L3 ANSWER 3 OF 8 CAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 2001:131163 CAPLUS

DOCUMENT NUMBER: 134:168379

TITLE: Preparation of time-specific controlled-release capsule formulations containing a swellable polymeric coating layers

INVENTOR(S): Busetti, Cesare; Crimella, Tiziano

PATENT ASSIGNEE(S): Italy

SOURCE: U.S., 11 pp., Cont.-in-part of U.S. 5,891,474.

CODEN: USXXAM

DOCUMENT TYPE: Patent

LANGUAGE: English

FAMILY ACC. NUM. COUNT: 2

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
US 6190692	B1	20010220	US 1997-991814	19971216
US 5891474	A	19990406	US 1997-790530	19970129
PRIORITY APPLN. INFO.: REFERENCE COUNT:	64	THERE ARE 64 CITED REFERENCES AVAILABLE FOR THIS RECORD. ALL CITATIONS AVAILABLE IN THE RE FORMAT		

AB The time-specific controlled-release capsule formulations comprise (a) a core contg. a liq. form of a pharmaceutically active agent to be delivered, and (b) a swellable polymeric coating layer substantially surrounding the core. The swellable polymeric coating layer delays the release of the pharmaceutically active agent from the core for a predetd. period of time dependent upon the thickness of the swellable polymeric coating layer. The swellable polymeric coating layer surrounding the core is provided by a new method which includes alternately (i) wetting the core with a binder soln., and (ii) coating the core with powd. polymeric particles a sufficient no. of times to produce a time-specific dosage formulation having the desired thickness of swellable polymeric coating layer. For example, 40 mg of verapamil HCl, 129 mg of dibasic calcium phosphate dihydrate, 20 mg of microcryst. cellulose, and 10 mg of sodium starch glycolate, were mixed thoroughly. Magnesium stearate (1 mg) is added and thoroughly mixed for another 5 min. The granular mixt. is formed into tablet cores of 6.8 mm diam., weighing 200 mg each using a rotary tablet press. The cores show a disintegration time lower than 5 min. in water, a Schleuninger hardness higher than 10 kp and a friability lower than 0.1 %. The cores are heated to 400.degree. and the coating layer is applied onto the cores in a two-step procedure, using an automatic coating pan. In the first step, the cores are wetted with a binder soln. contg. 5% Methocel E5, 10% polyvinylpyrrolidone, and 85% purified water. In the second step, the wetted cores were treated with a dry mixt. including 90% Methocel K15M, 9% talc and 1% colloidal silica. Steps 1 and 2 are repeated until a wt. gain corresponding to 50% of total tablet wt. is achieved. The coated tablets showed a dissoln. time lag in excess of 300 min., followed by a quick disintegration of the tablet.

IT Biopolymers

Gelatins, biological studies

Nucleic acids

Peptides, biological studies

Phosphatidylcholines, biological studies

Polymers, biological studies

Polyoxyalkylenes, biological studies

Steroids, biological studies

RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)
(prepn. of time-specific controlled-release capsules comprising
drug-contg. core and swellable polymeric coatings)

IT 7631-86-9, **Silica**, biological studies

RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)
(colloidal; prepn. of time-specific controlled-release capsules
comprising drug-contg. core and swellable polymeric coatings)

IT 50-70-4, Sorbitol, biological studies 56-81-5, Glycerol, biological
studies 57-88-5, Cholesterol, biological studies 58-95-7, Tocopherol
acetate 63-42-3, Lactose 64-17-5, Ethyl alcohol, biological studies
69-89-6D, **Xanthine**, derivs. 79-10-7D, Acrylic acid, esters,
polymers 79-41-4D, Methacrylic acid, esters, copolymers 110-16-7,
Maleic acid, biological studies 110-27-0, Isopropyl myristate 111-90-0
152-11-4, Verapamil hydrochloride 361-09-1, Sodium cholate 557-04-0,
Magnesium stearate 846-49-1, Lorazepam 7789-77-7, Calcium phosphate
dihydrate 9000-01-5, Arabic gum 9000-30-0, Guar gum 9000-69-5,
Pectin 9002-89-5, Polyvinyl alcohol 9003-01-4, Poly(acrylic acid)
9003-39-8, Polyvinylpyrrolidone 9004-10-8, Insulin, biological studies
9004-62-0, Hydroxyethyl cellulose 9004-64-2, Hydroxypropyl cellulose
9004-65-3, Hydroxypropyl methyl cellulose 9063-38-1, Sodium starch
glycolate 9087-70-1, Aprotinin 11138-66-2, Xanthan gum 14807-96-6,
Talcum, biological studies 15307-79-6, Diclofenac sodium 16051-77-7,
Isosorbide-5-mononitrate 18641-57-1, Glyceryl behenate 22260-51-1,
Bromocryptine mesylate 25322-68-3, PEG 33419-42-0, Etoposide
47931-85-1, Salmon calcitonin 59865-13-3, Cyclosporin A 65381-09-1,
Caprylic/capric triglyceride 65381-09-1D, Caprylic/capric triglyceride,
ethoxylated 83138-62-9, Polyglyceryl isostearate 106392-12-5,
Poloxamer

RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)
(prepn. of time-specific controlled-release capsules comprising
drug-contg. core and swellable polymeric coatings)

L3 ANSWER 4 OF 8 CAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 2000:384565 CAPLUS

DOCUMENT NUMBER: 133:28236

TITLE: Methods and compositions for performing an array of
chemical reactions on a support surface

INVENTOR(S): Zebala, John A.

PATENT ASSIGNEE(S): Syntrix Biochip, Inc., USA

SOURCE: PCT Int. Appl., 157 pp.

CODEN: PIXXD2

DOCUMENT TYPE: Patent

LANGUAGE: English

FAMILY ACC. NUM. COUNT: 4

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
WO 2000033084	A2	20000608	WO 1999-US28021	19991123
WO 2000033084	A3	20000810		
W:	AE, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, CA, CH, CN, CR, CU, CZ, DE, DK, DM, EE, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MA, MD, MG, MK, MN, MW, MX, NO, NZ, PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TR, TT, TZ, UA, UG, US, UZ, VN, YU, ZA, ZW, AM, AZ, BY, KG, KZ, MD, RU, TJ, TM			
RW:	GH, GM, KE, LS, MW, SD, SL, SZ, TZ, UG, ZW, AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, BF, BJ, CF, CG, CI, CM, GA, GN, GW, ML, MR, NE, SN, TD, TG			
AU 2000018317	A5	20000619	AU 2000-18317	19991123
EP 1163374	A2	20011219	EP 1999-961813	19991123
R:	AT, BE, CH, DE, DK, ES, FR, GB, GR, IT, LI, LU, NL, SE, MC, PT, IE, SI, LT, LV, FI, RO			

JP 2002531470	T2	20020924	JP 2000-585669	19991123
PRIORITY APPLN. INFO.:			US 1998-110527P	P 19981201
			US 1999-326479	A 19990604
			WO 1999-US28021	W 19991123

ST support array chem reaction photoresist; ligand array; DNA hybridization
immobilized probe; ACE inhibitor screening enalaprilat analog solid phase
synthesis; **nucleic acid** array

IT **Nucleic acid** hybridization
(DNA-DNA; methods and compns. for performing arrays of chem. reactions
on support surfaces using photoresists)

IT Adhesives

Analysis

Chromatography

DNA sequence analysis

Diagnosis

Drug screening

Electrophoresis

Human immunodeficiency virus

Indicators

Mass spectrometry

NMR spectroscopy

Negative photoresists

Nucleic acid hybridization

PCR (polymerase chain reaction)

Photoresists

Positive photoresists

Protein sequence analysis

RNA sequence analysis

Radiation

Reactors

Solvents

Surface

Synthesis
(methods and compns. for performing arrays of chem. reactions on
support surfaces using photoresists)

IT Probes (**nucleic acid**)
RL: ARG (Analytical reagent use); ANST (Analytical study); USES (Uses)
(methods and compns. for performing arrays of chem. reactions on
support surfaces using photoresists)

IT Peptide **nucleic acids**
RL: ARG (Analytical reagent use); DEV (Device component use); PEP
(Physical, engineering or chemical process); RCT (Reactant); SPN
(Synthetic preparation); ANST (Analytical study); PREP (Preparation); PROC
(Process); RACT (Reactant or reagent); USES (Uses)
(methods and compns. for performing arrays of chem. reactions on
support surfaces using photoresists)

IT **Nucleic acids**
Polynucleotides

Proteins, general, reactions

Reagents
RL: ARG (Analytical reagent use); DEV (Device component use); RCT
(Reactant); ANST (Analytical study); RACT (Reactant or reagent); USES
(Uses)
(methods and compns. for performing arrays of chem. reactions on
support surfaces using photoresists)

IT Peptides, reactions
RL: ARG (Analytical reagent use); DEV (Device component use); PEP
(Physical, engineering or chemical process); RCT (Reactant); ANST
(Analytical study); PROC (Process); RACT (Reactant or reagent); USES
(Uses)
(**nucleic acid** mimics; methods and compns. for
performing arrays of chem. reactions on support surfaces using
photoresists)

IT 51-20-7, 5-Bromouracil 51-21-8, 5-Fluorouracil 58-63-9, Inosine

65-71-4, Thymine 66-22-8, Uracil, uses 66-22-8D, Uracil, pseudo-,
 derivs., uses 68-94-0, Hypoxanthine 69-89-6, **Xanthine**
 71-30-7, Cytosine 73-24-5, Adenine, uses 73-40-5, Guanine 141-90-2,
 Thiouracil 333-49-3, 2-Thiocytosine 443-72-1 504-07-4, Dihydrouracil
 554-01-8, 5-Methylcytosine 578-76-7, 7-Methylguanine 591-28-6,
 4-Thiouracil 636-26-0, 5-Methyl-2-thiouracil 696-07-1, 5-Iodouracil
 938-85-2, 1-Methylguanine 1445-08-5, 2-Methyladenine 1445-15-4
 1500-85-2, 7-Deazaadenine 1820-81-1, 5-Chlorouracil 1904-98-9,
 2,6-Diaminopurine 2140-73-0, 1-Methylinosine 2365-40-4,
 N6-Isopentenyladenine 4776-08-3, 3-Methylcytosine 5142-22-3,
 1-Methyladenine 6623-81-0, 5-Methoxyuracil 7355-55-7, 7-Deazaguanine - - - - -
 10030-78-1 14631-20-0 14886-75-0 20758-33-2 31458-37-4
 72704-66-6 273752-46-8 273752-47-9 273752-48-0 273752-50-4
 273752-52-6

RL: DEV (Device component use); PRP (Properties); USES (Uses)
 (array of nucleobase polymers contg.; methods and compns. for
 performing arrays of chem. reactions on support surfaces using
 photoresists)

IT 7631-86-9, **Silica**, uses

RL: DEV (Device component use); USES (Uses)
 (methods and compns. for performing arrays of chem. reactions on
 support surfaces using photoresists)

L3 ANSWER 5 OF 8 CAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 1995:507908 CAPLUS

DOCUMENT NUMBER: 122:265933

TITLE: Preparation of pyranose nucleoside derivatives as
 antiviral and antitumor agents

INVENTOR(S): Waga, Toshiaki; Meguro, Hiromu; Oorui, Hiroshi

PATENT ASSIGNEE(S): Asahi Breweries Ltd, Japan

SOURCE: Jpn. Kokai Tokkyo Koho, 10 pp.

CODEN: JKXXAF

DOCUMENT TYPE: Patent

LANGUAGE: Japanese

FAMILY ACC. NUM. COUNT: 1

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
JP 06263793	A2	19940920	JP 1993-139791	19930304
PRIORITY APPLN. INFO.:			JP 1993-139791	19930304

OTHER SOURCE(S): MARPAT 122:265933

AB The title compds. (I; B = adenine, guanine, thymine, uracil, cytosine,
 hypoxanthine, **xanthine**, 5-methylcytosine, 4-ethoxy-5-methyl-2-
 oxopyrimidine, 4-isopropoxy-5-methyl-2-oxopyrimidine, 5-methyl-2-
 oxopyrimidine; R1, R2 = H, OH; or R1R2 = bond; R3 = Q wherein n = 0,1,3;
 R4 = H, lower alkoxy) or pharmacol. acceptable esters, ethers, or salts
 thereof are prepd. as antiviral and antitumor agents, particularly
 potential anti-HIV agents (no data), are prepd. Thus, 2.0 g adenine and
 2.0 g K2CO3 were suspended in 100 mL DMF and after stirring at 80.degree.
 for 1 h, 2.0 g 18-crown-6 ether and Me 2,3-anhydro-4,6-O-benzylidene-
 .alpha.-D-allopyranoside were added followed by stirring the resulting
 mixt. at 120.degree. for 16 h to give, after **silica** gel
 chromatog., 91% adenylaltropyranoside deriv. (II; RR = CHPh).

IT 3150-15-0

RL: RCT (Reactant); RACT (Reactant or reagent)
 (condensation with **nucleic acid** bases in prepn. of
 pyranose nucleoside derivs. as antiviral and antitumor agents)

L3 ANSWER 6 OF 8 CAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 1993:423983 CAPLUS

DOCUMENT NUMBER: 119:23983

TITLE: Optimized separation of purine bases and nucleosides
 in human cord plasma by capillary zone electrophoresis

AUTHOR(S): Grune, Tilman; Ross, Gordon A.; Schmidt, Heike; Siems, Werner; Perrett, David
CORPORATE SOURCE: Institute of Biochemistry, Medical Faculty (Charite), Humboldt University, Hessische Strasse 3-4, Berlin, O-1040, Germany
SOURCE: Journal of Chromatography (1993), 636(1), 105-11
CODEN: JOCRAM; ISSN: 0021-9673
DOCUMENT TYPE: Journal
LANGUAGE: English

AB An optimized sepn. of the main purine compds. of human serum by capillary zone electrophoresis is presented. Sepns. were performed in an uncoated ~~silica~~ capillary (44 cm .times. 75 .mu.m inner diam., 37 cm to window) on a SpectraPhoresis 1000 system with UV detection. The sepn. of adenine (Ade), adenosine (Ado), guanine (Gua), guanosine (Guo), hypoxanthine (Hyp), inosine (Ino), **xanthine** (Xan), and uric acid (UA) was optimized with respect to pH, temp., applied potential, and hydrodynamic injection time. Optimum conditions were 20 mM borate buffer (pH 9.4), 37.degree., 20 kV and 9 s load and detection at 260 nm. Linearity extended from 1 to 125 .mu.M. The sensitivity of the method was 0.5 .mu.M, which is adequate for measuring Ade, Gua, Hyp, and UA in plasma samples. Plasma samples from newborns were pptd. with an equal vol. of HClO4 (7%, vol./vol.), the supernatant was adjusted to neutral pH with K carbonate and, before injection, the sample was alkalized with NaOH. The method presented here allows the detn. of Ade, Guo, Hyp, and UA. The levels of the detd. purines were compared in samples from control newborns, preterm babies, and newborns with asphyxia or acidic serum pH values.

IT **Nucleic acid** bases
Nucleosides, analysis
RL: PROC (Process)

(purine, sepn. of, of newborn plasma by capillary zone electrophoresis)
IT 58-61-7, Adenosine, analysis 58-63-9, Inosine 68-94-0, Hypoxanthine 69-89-6, **Xanthine** 69-93-2, Uric acid, analysis 73-24-5, Adenine, analysis 73-40-5, Guanine 118-00-3, Guanosine, analysis
RL: ANT (Analyte); ANST (Analytical study)
(detn. and sepn. of, from purine bases and nucleosides of newborn plasma by capillary electrophoresis)

L3 ANSWER 7 OF 8 CAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 1984:628477 CAPLUS
DOCUMENT NUMBER: 101:228477
TITLE: Separation of nucleobases on polar amino cyano high-performance liquid chromatography columns
AUTHOR(S): Joshua, Henry; Goetz, Michael
CORPORATE SOURCE: Merck Sharp and Dohme Res. Lab., Rahway, NJ, 07065, USA
SOURCE: Journal of Chromatography (1984), 303(1), 185-9
CODEN: JOCRAM; ISSN: 0021-9673
DOCUMENT TYPE: Journal
LANGUAGE: English

AB The use of a polar amino cyano (PAC) HPLC column in conjunction with buffered water-acetonitrile (8:92) eluents affords a viable alternative to the std. reversed-phase, **silica** gel and ion-exchange methods for the chromatog. sepn. of nucleobases. A Whatman PAC column (Partisil PXS 5/25 PAC, 25 cm .times. 4.6 mm, 5 .mu.m particle size) was used equipped with an Upchurch C-130 precolumn packed with Whatman Co-Pell PAC 30-38 .mu.m particles. Seven nucleobases (adenine [73-24-5], cytosine [71-30-7], guanine [73-40-5], hypoxanthine [68-94-0], thymine [65-71-4], uracil [66-22-8], and **xanthine** [69-89-6]) were successfully sepd. Changes in the eluent pH values were found to affect the selectivity and capacity factors for the nucleobases. Thus effects were esp. pronounced with **xanthine**. The method is useful for the sepn. of nucleobases from fermn. broth exts.

ST **nucleic acid** base sepn polar HPLC; amino cyano

chromatog adenine cytosine guanine
IT Fermentation
 (nucleic acid bases, sepn. after, by polar
 aminocyano high-performance liq. chromatog.)
IT **Nucleic acids**
 RL: PROC (Process)
 (bases, sepn. of, from fermn. broth by polar amino cyano
 high-performance liq. chromatog.)
IT Chromatography, column and liquid
 (high-performance, **nucleic acid** bases sepn. from
 fermn. broth by, on polar amino cyano stationary phase) - - - - -

L3 ANSWER 8 OF 8 CAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 1968:48556 CAPLUS

DOCUMENT NUMBER: 68:48556

TITLE: Mode of action of the chemosterilants,
2-imidazolidinone and 4-imidazolin-2-one, in the
housefly and in the large milkweed bug

AUTHOR(S): Schaefer, Charles Herbert

CORPORATE SOURCE: Shell Develop. Co., Modesto, CA, USA

SOURCE: Life Sciences (1967), 6(24), 2677-83

CODEN: LIFSAK; ISSN: 0024-3205

DOCUMENT TYPE: Journal

LANGUAGE: English

AB An anal. method was developed for detecting 2-imidazolidinone (I) and 4-imidazolin-2-one (II) in insect tissues and feces after injection of 5 .mu.g. I or II into female houseflies or large milkweed bugs, *Oncopeltus fasciatus*. Insects were homogenized and the homogenate or feces were extd. with MeOH. After evapn., the residue was added to a H2O-CHCl3 mixt., centrifuged, the aq. phase extd. with CHCl3, sepd., and evapd. onto silica gel. Exts. were subjected to column chromatog. and thin-layer chromatog. Excretion plus metabolism eliminated 84% of a 5-.mu.g. dose of I and 76% of II within 24 hrs. after injection into houseflies; these 2 compds. are temporary sterilants in this species. The nature of the metabolites of either compd. was not detd. Large milkweed bugs were apparently unable to detoxify either compd. and there was no trace of either in the feces at 48 hrs. after injection; both I and II produce permanent sterility in this species. An attempt was made to det. the mode of action of II in houseflies by feeding the lowest level that inhibited reproduction in the presence of potential reversers; none of the natural biochem. intermediates tested (vitamins, glycine, histamine, inosinic acid, orotic acid, cytidine, adenine, cytidylic acid, adenylic acid, deoxyadenosine, **xanthine**, RNA, oleic acid, .beta.-sitosterol, and cholesterol) had any effect on the sterilant activity of II. No synergism was apparent when either compd. was fed in the presence of 0.2% sesamex or when II was given in the presence of 1% H3BO3. II may inhibit the formation of complex mols. such as proteins or **nucleic acids** rather than that of simple mols.; the mode of action of II is complex and may also involve interference with endocrine regulation.

=> d his

(FILE 'HOME' ENTERED AT 16:47:21 ON 13 JUL 2003)

FILE 'BIOSIS, CAPLUS, BIOTECHNO' ENTERED AT 16:47:30 ON 13 JUL 2003

L1 437637 S SILICA

L2 1076 S L1 AND NUCLEIC ACID

L3 8 S L2 AND XANTHINE

=> s l1(w)gel

L4 88477 L1(W) GEL

=> s l4 and nucleic acid

2 FILES SEARCHED...

L5 306 L4 AND NUCLEIC ACID

=> s l4 and nucleic acid purification

L6 17 L4 AND NUCLEIC ACID PURIFICATION

=> d 1-17 l6 ibib kwic

L6 ANSWER 1 OF 17 BIOSIS COPYRIGHT 2003 BIOLOGICAL ABSTRACTS INC.

ACCESSION NUMBER: 2002:125720 BIOSIS

DOCUMENT NUMBER: PREV200200125720

TITLE: **Nucleic acid purification**
using **silica gel** and glass particles.

AUTHOR(S): Padhye, V. V; York, C.; Burkiewicz, A.

CORPORATE SOURCE: Madison, Wis. USA

ASSIGNEE: PROMEGA CORPORATION

PATENT INFORMATION: US 5808041 Sept. 15, 1998

SOURCE: Official Gazette of the United States Patent and Trademark
Office Patents, (Sept. 15, 1998) Vol. 1214, No. 3, pp.
3017.

ISSN: 0098-1133.

DOCUMENT TYPE: Patent

LANGUAGE: English

TI **Nucleic acid purification** using
silica gel and glass particles.

IT Miscellaneous Descriptors

BIOTECHNOLOGY; GLASS PARTICLE; **NUCLEIC ACID**
PURIFICATION; SILICA GEL

L6 ANSWER 2 OF 17 BIOSIS COPYRIGHT 2003 BIOLOGICAL ABSTRACTS INC.

ACCESSION NUMBER: 2002:82283 BIOSIS

DOCUMENT NUMBER: PREV200200082283

TITLE: **Nucleic acid purification** on
silica gel and glass mixtures.

AUTHOR(S): Padhye, V. V; York, C.; Burkiewicz, A.

CORPORATE SOURCE: Madison, Wis. USA

ASSIGNEE: PROMEGA CORPORATION

PATENT INFORMATION: US 5658548 Aug. 19, 1997

SOURCE: Official Gazette of the United States Patent and Trademark
Office Patents, (Aug. 19, 1997) Vol. 1201, No. 3, pp. 2038.
ISSN: 0098-1133.

DOCUMENT TYPE: Patent

LANGUAGE: English

TI **Nucleic acid purification** on **silica**
gel and glass mixtures.

IT Miscellaneous Descriptors

BIOTECHNOLOGY; GLASS; METHODS; **NUCLEIC ACID**
PURIFICATION; SILICA GEL

L6 ANSWER 3 OF 17 BIOSIS COPYRIGHT 2003 BIOLOGICAL ABSTRACTS INC.

ACCESSION NUMBER: 2001:371571 BIOSIS

DOCUMENT NUMBER: PREV200100371571

TITLE: Endotoxin reduction in **nucleic acid**
purification.

AUTHOR(S): Smith, Craig E.; Creswell, Donald A. (1); Bitner, Rex M.;
White, Douglas H.; Butler, Braeden L.; Lesley, Scott A.

CORPORATE SOURCE: (1) Cottage Grove, WI USA

ASSIGNEE: Promega Corporation

PATENT INFORMATION: US 6194562 February 27, 2001

SOURCE: Official Gazette of the United States Patent and Trademark
Office Patents, (Feb. 27, 2001) Vol. 1243, No. 4, pp. No
Pagination. e-file.
ISSN: 0098-1133.

DOCUMENT TYPE: Patent
LANGUAGE: English

TI Endotoxin reduction in **nucleic acid**
purification.

AB. . . solutions contaminated with endotoxins from external sources. The present method removes endotoxins from such solutions using silica-based materials, such as **silica gel** particles, magnetic silica particles, or diatomaceous earth. In a preferred aspect of the method of this invention, magnetic silica particles. . .

IT Major Concepts

Molecular Genetics (Biochemistry and Molecular Biophysics); Methods and Techniques

IT Chemicals & Biochemicals

endotoxins: removal; **nucleic acids:**
purification; plasmid DNA: isolation; silica-based materials

IT Methods & Equipment

endotoxin removal: purification method; **nucleic acid**
purification: purification method

L6 ANSWER 4 OF 17 CAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 2002:429115 CAPLUS

DOCUMENT NUMBER: 137:2749

TITLE: Purification of DNA sequencing reactions using silica magnetic particles

INVENTOR(S): Bjerke, Michael P.; Otto, Paul E.

PATENT ASSIGNEE(S): Promega Corporation, USA

SOURCE: PCT Int. Appl., 34 pp.

CODEN: PIXXD2

DOCUMENT TYPE: Patent

LANGUAGE: English

FAMILY ACC. NUM. COUNT: 1

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
WO 2002044414	A2	20020606	WO 2001-US43364	20011121
WO 2002044414	A3	20030522		

W: AE, AG, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, BZ, CA, CH, CN, CO, CR, CU, CZ, DE, DK, DM, DZ, EC, EE, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MA, MD, MG, MK, MN, MW, MX, MZ, NO, NZ, PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TR, TT, TZ, UA, UG, UZ, VN, YU, ZA, ZW, AM, AZ, BY, KG, KZ, MD, RU, TJ, TM

RW: GH, GM, KE, LS, MW, MZ, SD, SL, SZ, TZ, UG, ZM, ZW, AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, TR, BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, ML, MR, NE, SN, TD, TG

AU 2002017784	A5	20020611	AU 2002-17784	20011121
---------------	----	----------	---------------	----------

PRIORITY APPLN. INFO.:	US 2000-724169	A	20001128
	WO 2001-US43364	W	20011121

IT Primers (**nucleic acid**)

RL: REM (Removal or disposal); PROC (Process)
(**purifn.** of DNA sequencing reactions using silica magnetic particles)

IT **Silica gel**, uses

RL: DEV (Device component use); USES (Uses)
(silica magnetic particles of; **purifn.** of DNA sequencing reactions using silica magnetic particles)

L6 ANSWER 5 OF 17 CAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 2002:51620 CAPLUS

DOCUMENT NUMBER: 136:97266

TITLE: Isolating nucleic acids by selective adsorption and desorption onto silicon dioxide

INVENTOR(S): Weber, Martin; Singer, Thorsten; Cosaert, Sarah

PATENT ASSIGNEE(S): Qiagen G.m.b.H., Germany
 SOURCE: PCT Int. Appl., 21 pp.
 CODEN: PIXXD2
 DOCUMENT TYPE: Patent
 LANGUAGE: German
 FAMILY ACC. NUM. COUNT: 1
 PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
WO 2002004620	A2	20020117	WO 2001-EP8066	20010712
WO 2002004620	A3	20020718		
W: JP, US				
RW: AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, TR				
DE 10033991	A1	20020124	DE 2000-10033991	20000712
EP 1299531	A2	20030409	EP 2001-971766	20010712
R: AT, BE, CH, DE, DK, ES, FR, GB, GR, IT, LI, LU, NL, SE, MC, PT, IE, FI, CY, TR				

PRIORITY APPLN. INFO.: DE 2000-10033991 A 20000712
 WO 2001-EP8066 W 20010712

ST nucleic acid purifn silica sorbent; alkali
 halide alc nucleic acid purifn; chloride
 isopropanol ethanol nucleic acid purifn

IT Glass fibers, uses
 Silica gel, uses

RL: DEV (Device component use); USES (Uses)
 (as sorbent; isolating nucleic acids by selective adsorption and
 desorption onto silicon dioxide)

L6 ANSWER 6 OF 17 CAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 2000:824267 CAPLUS

DOCUMENT NUMBER: 133:360593

TITLE: pH-dependent ion exchange matrix and method of
 synthesis and use for isolation of nucleic acids
 INVENTOR(S): Smith, Craig E.; Holmes, Diana L.; Simpson, Daniel J.;
 Katzhendler, Jehoshua; Bitner, Rex M.; Grosch,
 Josephine C.

PATENT ASSIGNEE(S): Promega Corp., USA

SOURCE: PCT Int. Appl., 63 pp.
 CODEN: PIXXD2

DOCUMENT TYPE: Patent
 LANGUAGE: English

FAMILY ACC. NUM. COUNT: 1

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
WO 2000069872	A2	20001123	WO 2000-US12186	20000505
WO 2000069872	A3	20010215		
W: AE, AG, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, CA, CH, CN, CR, CU, CZ, DE, DK, DM, DZ, EE, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MA, MD, MG, MK, MN, MW, MX, NO, NZ, PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TR, TT, TZ, UA, UG, UZ, VN, YU, ZA, ZW, AM, AZ, BY, KG, KZ, MD, RU, TJ, TM				
RW: GH, GM, KE, LS, MW, SD, SL, SZ, TZ, UG, ZW, AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, BF, BJ, CF, CG, CI, CM, GA, GN, GW, ML, MR, NE, SN, TD, TG				
US 6310199	B1	20011030	US 1999-312172	19990514
EP 1179057	A2	20020213	EP 2000-935865	20000505
R: AT, BE, CH, DE, DK, ES, FR, GB, GR, IT, LI, LU, NL, SE, MC, PT, IE, SI, LT, LV, FI, RO				
US 2001014650	A1	20010816	US 2001-813077	20010320

PRIORITY APPLN. INFO.:

US 1999-312172 A 19990514
WO 2000-US12186 W 20000505

ST **nucleic acid purifn** pH dependent ion
exchange matrix; DNA RNA plasmid purifn pH dependent ion exchange matrix
IT Glass fibers, reactions
Silica gel, reactions
RL: DEV (Device component use); RCT (Reactant); RACT (Reactant or
reagent); USES (Uses)
(solid support; pH-dependent ion exchange matrix and method of
synthesis and use for isolation of nucleic acids)

L6 -- ANSWER 7 OF 17 -- CAPLUS -- COPYRIGHT 2003 ACS

ACCESSION NUMBER: 1999:691107 CAPLUS
DOCUMENT NUMBER: 131:296203
TITLE: Removal of endotoxins during purification of nucleic
acids from bacterial sources
INVENTOR(S): Smith, Craig E.; Creswell, Donald A.; Bitner, Rex M.;
White, Douglas H.; Butler, Braeden L.; Lesley, Scott
A.
PATENT ASSIGNEE(S): Promega Corporation, USA
SOURCE: PCT Int. Appl., 49 pp.
CODEN: PIXXD2
DOCUMENT TYPE: Patent
LANGUAGE: English
FAMILY ACC. NUM. COUNT: 2
PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
WO 9954340	A1	19991028	WO 1999-US8491	19990422
W:	AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, CA, CH, CN, CU, CZ, DE, DK, EE, ES, FI, GB, GE, GH, GM, HR, HU, ID, IL, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MD, MG, MK, MN, MW, MX, NO, NZ, PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TR, TT, UA, UG, UZ, VN, YU, ZW, AM, AZ, BY, KG, KZ, MD, RU, TJ, TM			
RW:	GH, GM, KE, LS, MW, SD, SL, SZ, UG, ZW, AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, BF, BJ, CF, CG, CI, CM, GA, GN, GW, ML, MR, NE, SN, TD, TG			
US 6194562	B1	20010227	US 1998-64449	19980422
CA 2329067	AA	19991028	CA 1999-2329067	19990422
AU 9936513	A1	19991108	AU 1999-36513	19990422
AU 740145	B2	20011101		
EP 1071695	A1	20010131	EP 1999-918650	19990422
R:	BE, CH, DE, FR, GB, IT, LI, NL			
JP 2002512252	T2	20020423	JP 2000-544678	19990422
US 6284470	B1	20010904	US 2000-645133	20000824

PRIORITY APPLN. INFO.:

US 1998-64449 A 19980422
WO 1999-US8491 W 19990422
US 1999-134156P P 19990514
US 1999-475958 A3 19991230REFERENCE COUNT: 5 THERE ARE 5 CITED REFERENCES AVAILABLE FOR THIS
RECORD. ALL CITATIONS AVAILABLE IN THE RE FORMAT

AB A novel method for removing contaminating endotoxins from nucleic acids, such as DNA, RNA, or hybrids during purifn. are described. Nucleic acid sources that can be used include, but are not limited to, lysates of Gram-neg. bacteria and nucleic acid solns. contaminated with endotoxins from external sources. The present method removes endotoxins from such solns. using silica-based materials, such as **silica gel** particles, magnetic silica particles, or diatomaceous earth as sorbents. Two sorbents are used, one specific for nucleic acids, and the other specific for the endotoxins in a two-stage process. The first step uses the endotoxin sorbent and the second uses the non-bound material from the first stage and the nucleic acid sorbent. Preferably, the nucleic acids are prepd. in endotoxin-free water without the use of chaotropic agents.

In a preferred aspect of the method of this invention, magnetic silica particles are used to isolate plasmid DNA from a lysate of Gram-neg. bacteria transformed with the plasmid DNA. The preferred sorbents are com. silica matrixes. Application of the disclosed method produces nucleic acids which are sufficiently free of endotoxin contamination to be useful for a variety of different practical applications. Optimization expts. using cleared lysates from Escherichia coli JM109 are reported. Chaotropic agents (guanidine thiocyanate) were shown to inhibit endotoxin binding to the sorbent when the sorbent was pre-equilibrated with the agent, but the effects were less severe when the sorbent was suspended in water and the lysate contained the chaotropic agent. Comparison with prior art methods of endotoxin removal shows that the method of the invention lowers endotoxin concn. by .gtoreq.10-fold over older methods.

ST endotoxin removal **nucleic acid purifn** silica
sorbent; plasmid purifn endotoxin removal silica sorbent

IT Diatomite
Glass, uses
Silica gel, uses
RL: DEV (Device component use); USES (Uses)
(as sorbent; removal of endotoxins during purifn. of nucleic acids from bacterial sources)

IT Denaturants
(chaotropic, in **nucleic acid purifn.**;
removal of endotoxins during purifn. of nucleic acids from bacterial sources)

IT Chelating agents
(in **nucleic acid purifn.**; removal of
endotoxins during purifn. of nucleic acids from bacterial sources)

IT Alcohols, uses
RL: MOA (Modifier or additive use); USES (Uses)
(in **nucleic acid purifn.**; removal of
endotoxins during purifn. of nucleic acids from bacterial sources)

IT Salts, uses
RL: MOA (Modifier or additive use); USES (Uses)
(non-chaotropic, in **nucleic acid purifn.**;
removal of endotoxins during purifn. of nucleic acids from bacterial sources)

L6 ANSWER 8 OF 17 CAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 1999:607125 CAPLUS
DOCUMENT NUMBER: 131:201804
TITLE: Superparamagnetic adsorbents for purification of
nucleic acids by solid-phase extraction
INVENTOR(S): Schubert, Frank; Wambutt, Rolf
PATENT ASSIGNEE(S): AGOWA Gesellschaft fuer Molekularbiologische
Technologie m.b.H., Germany
SOURCE: Ger. Offen., 6 pp.
CODEN: GWXXBX
DOCUMENT TYPE: Patent
LANGUAGE: German
FAMILY ACC. NUM. COUNT: 1
PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
DE 19912799	A1	19990916	DE 1999-19912799	19990310
PRIORITY APPLN. INFO.:			DE 1998-19810709	19980312

AB Superparamagnetic **silica gel** adsorbents are prepd. for
isolation and purifn. of nucleic acids, DNA and RNA from PCR products and
body fluids by solid-phase extn. The paramagnetic substances, e.g.,
magnetite (Fe₃O₄) or ferrites, can be pptd. in the pores of the adsorbent
or incorporated into the adsorbent matrix by mixing FexOy
(2.ltoreq.x.ltoreq.3.5, 3.ltoreq.y.ltoreq.4.5) particles into an aq.
alkali silicate soln. followed by addn. of a C1-6-carboxylic acid (e.g.,

glacial acetic acid) for prepn. of the **silica gel**.
 The **silica gel** may be functionalized. In an example,
 $\text{FeCl}_2 \cdot 4\text{H}_2\text{O}$ and $\text{FeCl}_3 \cdot 6\text{H}_2\text{O}$ were dissolved in water and added dropwise to 2M
 NaOH. The pptd. iron oxide (av. size 80 nm) was sepd. by centrifugation,
 washed with water, then used in the synthesis of **silica**
gel from Na silicate. The resulting paramagnetic **silica**
gel was used for purifn. of PCR fragments from plasmid pKS.

ST superparamagnetic adsorbent **nucleic acid**
 purifn; solid phase extn **nucleic acid**
 purifn; DNA purifn solid phase extn superparamagnetic adsorbent;
 RNA purifn solid phase extn superparamagnetic adsorbent

IT **Silica gel**, uses
 RL: IMF (Industrial manufacture); NUU (Other use, unclassified); SPN
 (Synthetic preparation); PREP (Preparation); USES (Uses)
 (paramagnetic; superparamagnetic adsorbents for purifn. of nucleic
 acids by solid-phase extn.)

IT DNA
Nucleic acids
 Oligonucleotides
 RNA
 RL: PUR (Purification or recovery); PREP (Preparation)
 (purifn. of; superparamagnetic adsorbents for purifn. of
 nucleic acids by solid-phase extn.)

L6 ANSWER 9 OF 17 CAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 1998:604699 CAPLUS

DOCUMENT NUMBER: 129:227825

TITLE: **Nucleic acid purification**
 using **silica gel** and glass
 particles

INVENTOR(S): Padhye, Vikas V.; York, Chuck; Burkiewicz, Adam

PATENT ASSIGNEE(S): Promega Corporation, USA

SOURCE: U.S., 13 pp., Cont.-in-part of U.S. Ser. No. 115,504,
 abandoned.

CODEN: USXXAM

DOCUMENT TYPE: Patent

LANGUAGE: English

FAMILY ACC. NUM. COUNT: 3

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
US 5808041	A	19980915	US 1995-485429	19950607
CA 2170604	AA	19950309	CA 1994-2170604	19940830
PRIORITY APPLN. INFO.: US 1993-115504			B2	19930830
REFERENCE COUNT: 6 THERE ARE 6 CITED REFERENCES AVAILABLE FOR THIS RECORD. ALL CITATIONS AVAILABLE IN THE RE FORMAT				

TI **Nucleic acid purification** using
silica gel and glass particles

AB The present invention provides compns. and methods for isolating nucleic
 acids with lengths greater than about 50 bases, from cells, gels, solns.
 and other media, in which nucleic acids occur in vivo or in vitro. The
 compns. of the invention are mixts. of the silica materials **silica**
gel and glass particles, particularly glass microfibers; such
 mixts. combined with chaotropic salts, such as guanidinium chloride or
 guanidinium thiocyanate; and suspensions of such mixts. in aq. solns. of
 chaotropic salts. In the methods of the invention, an aq. soln.
 comprising nucleic acid is mixed with an aq. soln. of chaotropic salts and
 the resulting soln. is contacted with a mixt. of the silica materials,
 whereupon the nucleic acid in the soln. binds to the silica materials.
 The chaotropic salts and components, other than the nucleic acid adsorbed
 to the silica materials, from the aq. soln. treated by the method of the
 invention are washed from the silica materials. Finally, the nucleic acid
 can be obtained by elution from the silica materials. The methods provide

nucleic acid in water or buffer, such as TE buffer, free of contamination by any salt or macromol. that would interfere with further processing or anal.

ST **nucleic acid purifn silica gel glass**

IT Salts, uses

RL: NUU (Other use, unclassified); USES (Uses)
(Chaotropic; **nucleic acid purifn. using silica gel and glass particles**)

IT Buffers

Cell

Chelating agents

Gels

Purification

(**nucleic acid purifn. using silica gel and glass particles**)

IT Glass, uses

Silica gel, uses

RL: NUU (Other use, unclassified); USES (Uses)
(**nucleic acid purifn. using silica gel and glass particles**)

IT DNA

RL: PUR (Purification or recovery); PREP (Preparation)
(**nucleic acid purifn. using silica gel and glass particles**)

IT Nucleic acids

RL: PUR (Purification or recovery); PREP (Preparation)
(**nucleic acid purifn. using silica gel and glass particles**)

IT RNA

RL: PUR (Purification or recovery); PREP (Preparation)
(**nucleic acid purifn. using silica gel and glass particles**)

IT 7732-18-5, Water, analysis

RL: ANT (Analyte); ANST (Analytical study)
(**nucleic acid purifn. using silica gel and glass particles**)

IT 50-01-1, Guanidinium hydrochloride 60-00-4, Edta, uses 67-42-5, Ethyleneglycolbis(.beta.-aminoethylether)-N,N,N',N'-tetraacetic acid 77-92-9, uses 593-84-0, Guanidinium thiocyanate 650-51-1, Sodium trichloroacetate 7601-89-0, Sodium perchlorate 7681-82-5, Sodium iodide, uses

RL: NUU (Other use, unclassified); USES (Uses)
(**nucleic acid purifn. using silica gel and glass particles**)

L6 ANSWER 10 OF 17 CAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 1997:640679 CAPLUS

DOCUMENT NUMBER: 127:259756

TITLE: Process and device for isolating nucleic acids

INVENTOR(S): Lange, Hans

PATENT ASSIGNEE(S): Innova Gesellschaft zur Entwicklung und Vermarktung
Innovativer Produkte m.b, Germany; Lange, Hans

SOURCE: PCT Int. Appl., 50 pp.

CODEN: PIXXD2

DOCUMENT TYPE: Patent

LANGUAGE: German

FAMILY ACC. NUM. COUNT: 1

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
WO 9734908	A1	19970925	WO 1997-DE517	19970314
W: JP, US				

RW: AT, BE, CH, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE
 DE 19610354 C1 19971120 DE 1996-19610354 19960315
 EP 891369 A1 19990120 EP 1997-920510 19970314
 EP 891369 B1 20010816

R: AT, CH, DE, FR, GB, IT, LI, SE
 JP 2001503730 T2 20010321 JP 1997-533034 19970314
 AT 204285 E 20010915 AT 1997-920510 19970314
 US 6071395 A 20000606 US 1999-142958 19990125
 US 6232464 B1 20010515 US 1999-384936 19990827

PRIORITY APPLN. INFO.:

DE 1996-19610354 A 19960315
 WO 1997-DE517 W 19970314
 US 1999-142958 B3 19990125

AB A device is disclosed for purifying and concg. nucleic acids from biol. fluids and suspensions which contain nucleic acids prior to, e.g., their anal. by PCR. In the device, a reaction chamber contg. an adsorbent (e.g., **silica gel**, glass particles, glass fibers, or ion exchanger) is connected to a discharge chamber and the nucleic acids can, by means of an electrophoresis device, be transferred from the reaction chamber to the discharge chamber and concd. After purifn. of the nucleic acids, processes such as hybridization, amplification, and chemiluminescence detection can be performed.

ST **nucleic acid purifn** adsorption
 electrophoresis app

IT Immunoglobulins
 RL: NUU (Other use, unclassified); SPN (Synthetic preparation); PREP (Preparation); USES (Uses)
 (G, biotinylated; **nucleic acids purifn.** and concn. with adsorption/electrophoresis app.)

IT Recombination, genetic
 (amplification; **nucleic acids purifn.** and concn. with adsorption/electrophoresis app.)

IT Thermal cycling
 (app.; **nucleic acids purifn.** and concn. with adsorption/electrophoresis app.)

IT Proteins, specific or class
 RL: DEV (Device component use); NUU (Other use, unclassified); USES (Uses)
 (ligand-binding; **nucleic acids purifn.** and concn. with adsorption/electrophoresis app.)

IT Adsorbents
 Adsorption apparatus
 Blood
 Blood plasma
 Body fluid
 Chemiluminescence spectroscopy
 Computer application
 Computer program
 Electric conductors
 Electrodes
 Electrophoresis apparatus
 Ion exchangers
 Magnetic field
 Magnetic particles
 Magnets
 Membranes, nonbiological
 Microtiter plates
 Nucleic acid hybridization
 Photomultipliers
 Pipets
 Pumps
 (**nucleic acids purifn.** and concn. with adsorption/electrophoresis app.)

IT Antibodies
 Antigens
 Biopolymers

Glass, uses
Glass fibers, uses
Ligands
Metals, uses
Oligonucleotides
Peptide nucleic acids
Receptors

Silica gel, uses

RL: DEV (Device component use); NUU (Other use, unclassified); USES (Uses)
(**nucleic acids purifn.** and concn. with
adsorption/electrophoresis app.)

IT DNA

Nucleic acids

RL: PUR (Purification or recovery); PREP (Preparation)
(**nucleic acids purifn.** and concn. with
adsorption/electrophoresis app.)

IT Plastics, uses

RL: DEV (Device component use); NUU (Other use, unclassified); USES (Uses)
(thermoplastics; **nucleic acids purifn.**
and concn. with adsorption/electrophoresis app.)

IT 7439-89-6, Iron, uses 7440-22-4, Silver, uses 7782-42-5, Graphite,
uses 9012-36-6, Agarose

RL: DEV (Device component use); NUU (Other use, unclassified); USES (Uses)
(**nucleic acids purifn.** and concn. with
adsorption/electrophoresis app.)

IT 117710-36-8DP, IgG conjugates

RL: NUU (Other use, unclassified); SPN (Synthetic preparation); PREP
(Preparation); USES (Uses)
(**nucleic acids purifn.** and concn. with
adsorption/electrophoresis app.)

L6 ANSWER 11 OF 17 CAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 1997:574440 CAPLUS

DOCUMENT NUMBER: 127:201746

TITLE: **Nucleic acid purification**

on **silica gel** and glass mixtures

INVENTOR(S): Padhye, Vikas V.; York, Chuck; Burkiewicz, Adam

PATENT ASSIGNEE(S): Promega Corp., USA

SOURCE: U.S., 12 pp., Cont.-in-part of U.S. Ser. No. 115,504,
abandoned.

CODEN: USXXAM

DOCUMENT TYPE: Patent

LANGUAGE: English

FAMILY ACC. NUM. COUNT: 3

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
US 5658548	A	19970819	US 1995-476849	19950607
US 5658548	C1	20010724		
CA 2170604	AA	19950309	CA 1994-2170604	19940830

PRIORITY APPLN. INFO.: US 1993-115504 B2 19930830

TI **Nucleic acid purification on silica
gel and glass mixtures**

AB Nucleic acids with lengths greater than about 50 bases are isolated from
cells, gels, solns. and other media, in which nucleic acids occur in vivo
or in vitro, by using mixt. of **silica gel** and glass
microfibers combined with chaotropic salts such as guanidinium chloride or
guanidinium thiocyanate. An aq. soln. comprising nucleic acid is mixed
with an aq. soln. of chaotropic salts and the resulting soln. is contacted
with the above silica-based mixt. whereupon the nucleic acid in the soln.
binds to the silica materials. The chaotropic salts and components, other
than the nucleic acid adsorbed to the silica materials, are washed from
the silica materials and the nucleic acid is obtained by elution. The

methods provide nucleic acid in water or buffer free of contamination by any salt or macromol. that would interfere with further processing or anal.

ST RNA isolation siliceous carrier binding chaotropic; glass RNA isolation chaotropic salt; **silica gel nucleic acid purifn**

IT Glass fibers, biological studies

Silica gel, biological studies

RL: BSU (Biological study, unclassified); BIOL (Biological study)

(**nucleic acid purifn. on silica gel and glass mixts.**)

IT DNA

RL: PUR (Purification or recovery); PREP (Preparation)

(**nucleic acid purifn. on silica gel and glass mixts.**)

IT RNA

RL: PUR (Purification or recovery); PREP (Preparation)

(**nucleic acid purifn. on silica gel and glass mixts.**)

IT 50-01-1, Guanidine hydrochloride 60-00-4, Edta, properties 67-42-5, Egta 139-33-3, Edta disodium salt 593-84-0, Guanidinium thiocyanate 9003-98-9, DNase 142298-75-7, RNase inhibitor

RL: PRP (Properties); TEM (Technical or engineered material use); USES (Uses)

(**nucleic acid purifn. on silica gel and glass mixts.**)

L6 ANSWER 12 OF 17 CAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 1997:178962 CAPLUS

DOCUMENT NUMBER: 126:168833

TITLE: Purification, stabilization, or isolation of nucleic acids from biological materials

INVENTOR(S): Mueller, Oliver; Deuter, Rainer

PATENT ASSIGNEE(S): Max-Planck-Gesellschaft Zur Foerderung Der Wissenschaften E.V., Germany

SOURCE: Ger. Offen., 6 pp.

CODEN: GWXXBX

DOCUMENT TYPE: Patent

LANGUAGE: German

FAMILY ACC. NUM. COUNT: 1

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
DE 19530132	A1	19970220	DE 1995-19530132	19950816
DE 19530132	C2	19980716		
CA 2228769	AA	19970227	CA 1996-2228769	19960814
WO 9707239	A1	19970227	WO 1996-EP3595	19960814
W: AU, BR, CA, JP, MX, US				
RW: AT, BE, CH, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE				
AU 9668216	A1	19970312	AU 1996-68216	19960814
AU 712331	B2	19991104		
EP 851937	A1	19980708	EP 1996-928466	19960814
EP 851937	B1	20020403		
R: AT, BE, CH, DE, DK, FR, GB, IT, LI, LU, NL, SE				
JP 11511020	T2	19990928	JP 1997-508945	19960814
AT 215611	E	20020415	AT 1996-928466	19960814
US 6084091	A	20000704	US 1998-11567	19980211

PRIORITY APPLN. INFO.: DE 1995-19530132 A 19950816

WO 1996-EP3595 W 19960814

AB The invention concerns the purifn., stabilization, and/or isolation of nucleic acids from, e.g., tissues, body fluids, plants, microorganisms, feces as well as foods, sewage sludge, wastewater, etc., by adding a carbohydrate-based adsorption matrix to the nucleic acid-contg. sample in

an appropriate buffer to bind contaminants or impurities. The carbohydrate-based adsorbent can contain, e.g., starch, cellulose, potato flour, etc. The impurities in a nucleic acid-contg. sample can be, e.g., dehydr. products of Hbs and or bile acids or their salts. The sep'd. nucleic acids can be treated with enzymes for amplification and/or restriction cleavage reactions. The method may be used to isolate or detect nucleic acids from stool samples as a diagnostic test for tumors of the digestive tract, and esp. of the pancreas or intestine, and for bacterial or viral infections. Reagent kits are also disclosed for the purifn. and stabilization of nucleic acids of biol. materials, and the kits contain buffer, adsorption matrix for binding impurities, mineral carriers (e.g., metal oxides, silica gel, zeolites, etc.), and/or org. carriers (e.g., modified latex, synthetic polymers, or their mixts.), and other necessary solns. and accessories. An example is given of the anal. of DNA of human stool samples, comparing the capacities of bovine serum albumin, cellulose, potato starch, and potato flour as adsorption matrix, and potato flour was best.

- ST biol material **nucleic acid purifn** adsorbent;
feces DNA analysis adsorbent potato flour; tumor diagnosis feces nucleic acid detection; infection diagnosis feces nucleic acid detection; diagnosis nucleic acid detection adsorbent; digestive tract cancer diagnosis DNA detection
- IT Hemoglobins
RL: ARU (Analytical role, unclassified); ANST (Analytical study)
(dehydr. products; **nucleic acids purifn.**
and stabilization and isolation from biol. materials)
- IT Potato (Solanum tuberosum)
Potato (Solanum tuberosum)
(flour; **nucleic acids purifn.** and
stabilization and isolation from biol. materials)
- IT Nucleic acid amplification (method)
(ligase chain reaction; **nucleic acids purifn.** and stabilization and isolation from biol. materials)
- IT Digestive tract
(neoplasm; **nucleic acids purifn.** and
stabilization and isolation from biol. materials)
- IT Nucleic acid amplification (method)
(nucleic acid base-specific amplification; **nucleic acids purifn.** and stabilization and isolation from biol. materials)
- IT Adsorbents
Animal tissue
Bacteria (Eubacteria)
Biological materials
Body fluid
Bone marrow
Diagnosis
Feces
Filters
Food analysis
Fossils
Frits
Infection
Intestine, neoplasm
Latex
Membranes, nonbiological
Microorganism
Mutation
Neoplasm
PCR (polymerase chain reaction)
Pancreas, neoplasm
Particles
Plant analysis
Plant tissue

Purification
Soil analysis
Virus
Wastewater treatment
Wastewater treatment sludge
(**nucleic acids purifn.** and stabilization
and isolation from biol. materials)

IT DNA
Nucleic acids
RL: ANT (Analyte); PEP (Physical, engineering or chemical process); PUR
(Purification or recovery); THU (Therapeutic use); ANST (Analytical
study); BIOL (Biological study); PREP (Preparation); PROC (Process); USES
(Uses)
(**nucleic acids purifn.** and stabilization
and isolation from biol. materials)

IT Bile acids
Bile salts
Fibers
Glass, analysis
RL: ARU (Analytical role, unclassified); ANST (Analytical study)
(**nucleic acids purifn.** and stabilization
and isolation from biol. materials)

IT Oxides (inorganic), analysis
Polymers, analysis
Silica gel, analysis
Zeolites (synthetic), analysis
RL: ARU (Analytical role, unclassified); NUJ (Other use, unclassified);
ANST (Analytical study); USES (Uses)
(**nucleic acids purifn.** and stabilization
and isolation from biol. materials)

IT Albumins, analysis
RL: ARU (Analytical role, unclassified); NUJ (Other use, unclassified);
PEP (Physical, engineering or chemical process); ANST (Analytical study);
PROC (Process); USES (Uses)
(**nucleic acids purifn.** and stabilization
and isolation from biol. materials)

IT Carbohydrates, analysis
RL: ARU (Analytical role, unclassified); PEP (Physical, engineering or
chemical process); ANST (Analytical study); PROC (Process)
(**nucleic acids purifn.** and stabilization
and isolation from biol. materials)

IT Gene
RL: ANT (Analyte); THU (Therapeutic use); ANST (Analytical study); BIOL
(Biological study); USES (Uses)
(oncogene; **nucleic acids purifn.** and
stabilization and isolation from biol. materials)

IT Flours and Meals
Flours and Meals
(potato flour; **nucleic acids purifn.** and
stabilization and isolation from biol. materials)

IT Nucleic acid amplification (method)
(self-sustained sequence replication; **nucleic acids
purifn.** and stabilization and isolation from biol. materials)

IT Gene, animal
RL: ANT (Analyte); THU (Therapeutic use); ANST (Analytical study); BIOL
(Biological study); USES (Uses)
(tumor suppressor; **nucleic acids purifn.**
and stabilization and isolation from biol. materials)

IT 7732-18-5, Water, analysis
RL: AMX (Analytical matrix); ANST (Analytical study)
(**nucleic acids purifn.** and stabilization
and isolation from biol. materials)

IT 14808-60-7, Quartz, analysis
RL: ARU (Analytical role, unclassified); ANST (Analytical study)

(nucleic acids purifn. and stabilization
and isolation from biol. materials)
IT 9005-25-8, Potato starch, analysis
RL: ARU (Analytical role, unclassified); NUU (Other use, unclassified);
PEP (Physical, engineering or chemical process); ANST (Analytical study);
PROC (Process); USES (Uses)

(nucleic acids purifn. and stabilization
and isolation from biol. materials)
IT 9004-34-6, Cellulose, analysis
RL: ARU (Analytical role, unclassified); PEP (Physical, engineering or
chemical process); ANST (Analytical study); PROC (Process)
(nucleic acids purifn. and stabilization
and isolation from biol. materials)

L6 ANSWER 13 OF 17 CAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 1995:837638 CAPLUS
DOCUMENT NUMBER: 123:283735
TITLE: Chromatographic removal of endotoxins from
macromolecules manufactured by fermentation
INVENTOR(S): Colpan, Metin; Moritz, Peter; Schorr, Joachim
PATENT ASSIGNEE(S): Qiagen GmbH, Germany
SOURCE: PCT Int. Appl., 19 pp.
CODEN: PIXXD2
DOCUMENT TYPE: Patent
LANGUAGE: German
FAMILY ACC. NUM. COUNT: 3
PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
WO 9521179	A1	19950810	WO 1995-EP391	19950203
W: AU, CA, JP, US				
RW: AT, BE, CH, DE, DK, ES, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE				
DE 4432654	A1	19960321	DE 1994-4432654	19940914
DE 4432654	C2	19980326		
CA 2182388	AA	19950810	CA 1995-2182388	19950203
CA 2182397	AA	19950810	CA 1995-2182397	19950203
AU 9515777	A1	19950821	AU 1995-15777	19950203
AU 693511	B2	19980702		
WO 9608500	A1	19960321	WO 1995-EP392	19950203
W: AU, CA, JP, US				
RW: AT, BE, CH, DE, DK, ES, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE				
AU 9516647	A1	19960329	AU 1995-16647	19950203
EP 775150	A1	19970528	EP 1995-907641	19950203
EP 775150	B1	19990428		
R: AT, BE, CH, DE, DK, FR, GB, IE, IT, LI, LU, NL, SE				
EP 781291	A1	19970702	EP 1995-908259	19950203
R: CH, DE, FR, GB, LI				
JP 09508407	T2	19970826	JP 1995-520391	19950203
AT 179425	E	19990515	AT 1995-907641	19950203
US 5747663	A	19980505	US 1996-687522	19960930
US 6274371	B1	20010814	US 1997-809072	19970619
US 2003036175	A1	20030220	US 2002-254845	20020926
PRIORITY APPLN. INFO.:			DE 1994-4403692	A 19940207
			DE 1994-4422291	A 19940625
			DE 1994-4431125	A 19940901
			DE 1994-4432654	A 19940914
			WO 1995-EP391	W 19950203
			WO 1995-EP392	W 19950203
			US 1996-687522	A1 19960930
			US 1998-26613	B1 19980220
			US 1999-443091	B3 19991118
IT Nucleic acids				
RL: PUR (Purification or recovery); PREP (Preparation)				

(purifn. of; chromatog. removal of endotoxins from macromols.
manufd. by fermn.)

IT Kieselguhr

Silica gel, uses

RL: NUU (Other use, unclassified); USES (Uses)

(sorbent in prepn. endotoxin-free DNA; chromatog. removal of endotoxins
from macromols. manufd. by fermn.)

L6 ANSWER 14 OF 17 CAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 1995:341134 CAPLUS

DOCUMENT NUMBER: 122:101132

TITLE: Chromatographic purification and separation of nucleic
acid mixtures

INVENTOR(S): Feuser, Petra; Hermann, Ralf; Schorr, Joachim; Colpan,
Metin; Bastian, Helge

PATENT ASSIGNEE(S): Diagen Institut fuer Molekularbiologische Diagnostik
GmbH, Germany

SOURCE: Ger. Offen., 9 pp.

CODEN: GWXXBX

DOCUMENT TYPE: Patent

LANGUAGE: German

FAMILY ACC. NUM. COUNT: 1

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
DE 4321904	A1	19950112	DE 1993-4321904	19930701
CA 2142910	AA	19950112	CA 1994-2142910	19940624
CA 2142910	C	20020827		
WO 9501359	A1	19950112	WO 1994-EP2056	19940624
W: CA, JP, US				
RW: AT, BE, CH, DE, DK, ES, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE				
EP 658164	A1	19950621	EP 1994-922869	19940624
EP 658164	B1	20010404		
R: AT, BE, CH, DE, DK, ES, FR, GB, IE, IT, LI, NL, PT, SE				
JP 08501321	T2	19960213	JP 1994-503247	19940624
AT 200293	E	20010415	AT 1994-922869	19940624
ES 2155477	T3	20010516	ES 1994-922869	19940624
US 6383393	B1	20020507	US 1996-392882	19960315

PRIORITY APPLN. INFO.: DE 1993-4321904 A 19930701

WO 1994-EP2056 W 19940624

AB Nucleic acids are sepd. and purified from a nucleic acid mixt. by
adsorption from a high-ionic-strength aq. soln. contg. 1-50 vol.% Cl-5
aliph. alc., PEG, hydrophobic inorg. and/or org. polymer, and/or Cl₃CCO₂H
onto a porous or nonporous mineral carrier comprising a metal oxide,
silica gel, glass, or zeolite, washing the adsorbent,
and eluting with a soln. of lower ionic strength. Thus, a tissue sample
was homogenized in a soln. contg. 4-8M chaotropic salt (e.g.
guanidine-HCl, guanidine isothiocyanate, NaI), an org. solvent (e.g. PhOH,
CHCl₃, Et₂O), and detergent, digested with protease, mixed. with 0.5 vol.
95-100% aliph. alc. or PEG, and centrifuged, and the supernatant was
passed through an appropriate membrane or gel matrix which was washed with
an aq. soln. contg. 100 mM NaCl, 10 mM Tris-HCl (pH 7.5), and 30-80% alc.
or PEG to remove impurities. Nucleic acids were then eluted with 10 mM
Tris-HCl (pH 9.0) or distd. water for use in PCR.

ST **nucleic acid purifn** chromatog; adsorption

nucleic acid purifn

IT Ceramic materials and wares

Glass, oxide

Glass fibers, uses

Membranes

Oxides, uses

Silica gel, uses

Zeolites, uses

RL: ARG (Analytical reagent use); ANST (Analytical study); USES (Uses)
(chromatog. purifn. and sepn. of nucleic acid mixts.)

L6 ANSWER 15 OF 17 CAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 1993:444668 CAPLUS

DOCUMENT NUMBER: 119:44668

TITLE: chromatography-based apparatus and method for
isolation and purification of nucleic acids

INVENTOR(S): Colpan, Metin

PATENT ASSIGNEE(S): Diagen Institut fuer Moleukularbiologische Diagnostik
G.m.b.H., Germany

SOURCE: Ger. Offen., 25 pp.

CODEN: GWXXBX

DOCUMENT TYPE: Patent

LANGUAGE: German

FAMILY ACC. NUM. COUNT: 1

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
DE 4139664	A1	19930603	DE 1991-4139664	19911202
WO 9311218	A1	19930610	WO 1992-EP2774	19921201
W: JP, US				
RW: AT, BE, CH, DE, DK, ES, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE				
WO 9311221	A1	19930610	WO 1992-EP2775	19921201
W: JP, US				
RW: AT, BE, CH, DE, DK, ES, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE				
EP 616638	A1	19940928	EP 1992-924636	19921201
EP 616638	B1	19960410		
R: BE, CH, DE, DK, FR, GB, IT, LI, LU, NL, SE				
EP 616639	A1	19940928	EP 1992-924637	19921201
EP 616639	B1	19981104		
R: BE, CH, DE, DK, FR, GB, IT, LI, LU, NL, SE				
JP 07501223	T2	19950209	JP 1993-509824	19921201
JP 3115324	B2	20001204		
EP 875271	A2	19981104	EP 1998-107576	19921201
EP 875271	A3	20010425		
R: BE, CH, DE, DK, FR, GB, IT, LI, LU, NL, SE				
JP 2001095572	A2	20010410	JP 2000-223370	19921201
US 6277648	B1	20010821	US 1994-253152	19940602
US 2001047966	A1	20011206	US 2001-900199	20010709
PRIORITY APPLN. INFO.:			DE 1991-4139664	A 19911202
			EP 1992-924637	A3 19921201
			JP 1993-509824	A3 19921201
			WO 1992-EP2774	W 19921201
			WO 1992-EP2775	W 19921201
			US 1994-253152	A1 19940602

AB Nucleic acids are isolated from cells by disrupting the cells, removing cell debris preferably by filtration, adsorbing the nucleic acids on an anion exchanger in low-ionic-strength buffer, desorbing with a buffer of high ionic strength, adsorbing on an inorg. carrier at high ionic strength, and desorbing with water or low-ionic-strength buffer. A column for the process contains the 2 adsorbents in sep. segments which may be sepd. by a porous sintered glass or ceramic disk or a membrane. Thus, Escherichia coli HB 101 cells contg. plasmid pUC18 were lysed with NaOH-SDS, the lysate was centrifuged, and the supernatant was applied to a centrifugal extn. column contg. a DEAE anion exchanger and silica gel. After centrifugation, the anion exchanger segment was washed to remove RNA and proteins and eluted with 7M NaClO4-15% EtOH-10 mM NaOAc (pH 7.0) directly onto the silica gel segment. After washing the silica gel, DNA was eluted with 10 mM Tris-HCl-1 mM EDTA (pH 8.0). This DNA could be used directly for restriction cleavage, labeling, sequencing, or amplification.

ST nucleic acid purifn chromatog; DNA purifn

chromatog
IT Glass, oxide
Kaolin, uses
Silica gel, uses
Zeolites, uses
RL: ANST (Analytical study)
(nucleic acid purifn. by adsorption
chromatog. on, ionic strength in relation to)
IT Deoxyribonucleic acids
Nucleic acids
RL: PUR (Purification or recovery); PREP (Preparation)
(purifn. of, by chromatog. on ion exchanger and inorg.
adsorbent)
IT 1314-23-4, Zirconium oxide, properties 1344-28-1, Aluminum oxide,
properties 13463-67-7, Titanium dioxide, properties
RL: PRP (Properties)
(nucleic acid purifn. by adsorption
chromatog. on, ionic strength in relation to)
L6 ANSWER 16 OF 17 CAPLUS COPYRIGHT 2003 ACS
ACCESSION NUMBER: 1991:118112 CAPLUS
DOCUMENT NUMBER: 114:118112
TITLE: Method for purifying nucleic acids using an adsorbent
to remove contaminating proteins
INVENTOR(S): McCormick, Randy Miles
PATENT ASSIGNEE(S): du Pont de Nemours, E. I., and Co., USA
SOURCE: PCT Int. Appl., 30 pp.
CODEN: PIXXD2
DOCUMENT TYPE: Patent
LANGUAGE: English
FAMILY ACC. NUM. COUNT: 1
PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
WO 9010637	A1	19900920	WO 1989-US902	19890310
W: AU, DK, FI, JP, NO				
RW: AT, BE, CH, DE, FR, GB, IT, LU, NL, SE				
EP 462100	A1	19911227	EP 1989-903325	19890310
EP 462100	B1	19930526		
R: AT, BE, CH, DE, FR, GB, IT, LI, LU, NL, SE				
JP 04503942	T2	19920716	JP 1989-503163	19890310
JP 2703636	B2	19980126		
AU 8933525	A1	19921224	AU 1989-33525	19890310
AU 632284	B2	19920716		
AT 89827	E	19930615	AT 1989-903325	19890310
DK 9101561	A	19910905	DK 1991-1561	19910905
NO 9103545	A	19910909	NO 1991-3545	19910909
NO 177856	B	19950828		
NO 177856	C	19951206		
FI 95270	B	19950929	FI 1991-4240	19910909
FI 95270	C	19960110		

PRIORITY APPLN. INFO.: EP 1989-903325 19890310
WO 1989-US902 19890310

AB Proteins are sepd. from nucleic acids by contacting the mixt. with a solid
phase extn. material capable of binding proteins and then isolating the
unbound fraction contg. the nucleic acids. **Silica gel**
was treated with HF and rehydrated. The rehydrated gel removed all of the
PST, EcoRI, and Sall restriction enzymes from pBR322 DNA, enzyme, and
buffer mixts. before the enzymes could cut the DNA. DNA recovery from
various rehydrated **silica gels** was 76.1-92.9%.
ST **nucleic acid purifn** protein adsorption; DNA
purifn silica gel; restriction enzyme removal
silica gel

IT Bacteria
Escherichia coli
(DNA of, purifn. of, from blood, rehydrated **silica gel** removal of proteins in)

IT Blood
(bacterial DNA purifn. from, rehydrated **silica gel** removal of proteins in)

IT Adsorbents
(for protein removal in **nucleic acid purifn** .)

IT **Silica gel**, preparation
RL: SPN (Synthetic preparation); PREP (Preparation)
(prepn. of silane-treated, and use in protein removal from nucleic acid)

IT Deoxyribonucleic acids
Nucleic acids
RL: PUR (Purification or recovery); PREP (Preparation)
(purifn. of, adsorbents for protein removal in)

IT **Silica gel**, compounds
RL: ANST (Analytical study)
(rehydrated, for protein removal in **nucleic acid purifn.**)

IT Albumins, uses and miscellaneous
RL: REM (Removal or disposal); PROC (Process)
(removal of, from DNA, silane-treated **silica gel** for)

IT Plasmid and Episome
(pBR322, DNA of, restriction enzymes removal from, rehydrated **silica gel** in)

IT 7647-01-0, Hydrochloric acid, uses and miscellaneous 7664-39-3,
Hydrofluoric acid, uses and miscellaneous 7697-37-2, Nitric acid, uses
and miscellaneous
RL: ANST (Analytical study)
(in rehydrated **silica gel** prepn. for protein removal from DNA)

IT 9001-78-9, Alkaline phosphatase 80498-17-5, EcoRI Restriction enzyme
81295-32-1, PST restriction enzyme 81295-38-7, Sali Restriction enzyme
RL: REM (Removal or disposal); PROC (Process)
(removal of, from DNA, rehydrated **silica gel** for)

L6 ANSWER 17 OF 17 CAPLUS COPYRIGHT 2003 ACS

ACCESSION NUMBER: 1989:474361 CAPLUS

DOCUMENT NUMBER: 111:74361

TITLE: Anion exchanger on porous **silica gel**
matrix and its use in chromatographic purification of
long-chain nucleic acids

INVENTOR(S): Henco, Karsten; Stichel, Arndt; Colpan, Metin

PATENT ASSIGNEE(S): DIAGEN Institut fuer Molekularbiologische Diagnostik
G.m.b.H., Fed. Rep. Ger.

SOURCE: Ger. Offen., 9 pp.

CODEN: GWXXBX

DOCUMENT TYPE: Patent

LANGUAGE: German

FAMILY ACC. NUM. COUNT: 1

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
DE 3639949	A1	19880609	DE 1986-3639949	19861122
EP 268946	A2	19880601	EP 1987-116713	19871112
EP 268946	A3	19900314		
EP 268946	B1	19930915		
R: AT, BE, CH, DE, FR, GB, IT, LI, LU, NL, SE				
AT 94553	E	19931015	AT 1987-116713	19871112

CA 1339772	A1	19980324	CA 1987-552250	19871119
US 5057426	A	19911015	US 1987-123698	19871123
JP 63150294	A2	19880622	JP 1987-295947	19871124
JP 07013077	B4	19950215		

PRIORITY APPLN. INFO.:	DE 1986-3639949	19861122
	EP 1987-116713	19871112

TI Anion exchanger on porous **silica gel** matrix and its use in chromatographic purification of long-chain nucleic acids

AB Long-chain nucleic acids are purified from other substances in exts. of gently disrupted cells, body fluids, viruses, etc. by binding to a porous matrix of modified **silica gel** particles, washing out other substances, and elution of the nucleic acids. The matrix comprises silanized **silica gel** particles 15-250 .mu.m in size with pores 100-2500 nm in diam. bearing an anion exchanger, esp. N,N-dimethylaminoethanol linked to the matrix via .gamma.-glycidyloxypropyltrimethoxysilane. A suspension of lysed Escherichia coli cells contg. phage .lambda. was centrifuged and the supernatant was passed through a 0.45-.mu.m filter to remove intact cells and cell debris. The suspension was passed through a cartridge contg. the above modified **silica gel** to remove cellular DNA. The phage particles in the eluate were disrupted with 4M urea and passed through a cartridge contg. the above anion exchanger-modified **silica gel**, which bound the phage DNA. The cartridge was washed with 50 mM Tris-HCl buffer (pH 7.5) contg. 0.8M NaCl and 1 mM EDTA, and the phage .lambda. DNA was eluted with the same buffer contg. 1.2M NaCl and 1 mM EDTA and desalted by dialysis or pptn. with EtOH, PEG, or iso-PrOH.

ST **nucleic acid purifn silica gel**; DNA purifn anion exchanger **silica gel**

IT Liver, composition
Sperm
(DNA of, purifn. of, on porous **silica gel** anion exchanger)

IT Anion exchangers
(immobilized, on porous silanized **silica gel**, for **nucleic acid purifn.**)

IT Animal cell
Animal tissue
Bacteria
Body fluid
Plant cell
Plant tissue
Virus
(nucleic acids of, purifn. of, on porous **silica gel** anion exchanger)

IT Plasmid and Episome
(purifn. of, on porous **silica gel** anion exchanger)

IT Deoxyribonucleic acids
Nucleic acids
RL: PUR (Purification or recovery); PREP (Preparation)
(purifn. of, on porous **silica gel** anion exchanger)

IT **Silica gel**, compounds
RL: ANST (Analytical study)
(silanized, anion exchanger derivs., **nucleic acids purifn.** on)

IT Virus, animal
(cytomegalo-, DNA of, purifn. of, on porous **silica gel** anion exchanger)

IT Virus, bacterial
(lambda, DNA of, purifn. of, on porous **silica gel** anion exchanger)

IT Virus, animal
(papilloma, DNA of, purifn. of, on porous **silica gel** anion exchanger)